Influenza A/B Antigen Rapid Test (Nasal)

A rapid test for the qualitative detection of Influenza type A virus and Influenza type B virus nucleocapsid antigens in anterior nasal swab specimens.
For professional *in vitro* diagnostic use only.
Please read the instructions for use carefully before using.

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25 Tests/Kit

【INTENDED USE】

The GeneSign Influenza A/B Antigen Rapid Test is a lateral flow immunoassay intended for the *in vitro* rapid, simultaneous qualitative detection and differentiation of nucleocapsid antigen from Influenza type A virus and Influenza type B virus directly from anterior nasal swab specimens obtained from symptomatic patients who are suspected of respiratory Influenza type A or B viral infection by their healthcare providers. The test only provides preliminary analysis results for medical professionals' reference but not critical diagnosis criteria. Any reactive specimen with this test kit must be analyzed and confirmed with alternative testing method(s) and clinical findings. Final medical treatment and management of patients must be followed up with doctor's diagnosis and advice.

【SUMMARY】

Influenza is a highly contagious acute viral infection of the respiratory tract. It is a communicable disease easily transmitted from person to person through aerosol droplets excreted when sneezing and coughing. Common symptoms include high fever, chills, headache, cough, sore throat and malaise. The influenza A virus is more prevalent and is the primary pathogen associated with serious epidemics. The influenza B virus causes a disease that is generally not as severe as that caused by the influenza A virus.

An accurate diagnosis of influenza based on clinical symptoms is difficult because the initial symptoms of influenza are similar to many other illnesses. Therefore, it can be confirmed only by laboratory diagnostic testing. Early differential diagnosis of influenza type A or type B can allow for proper treatment with appropriate antiviral therapy while reducing the incidence of inappropriate treatment with antibiotics. Early diagnosis and treatment are of particular value in a clinical setting where an accurate diagnosis can assist the healthcare professionals with the management of influenza patients who are at risk for complications.

The GeneSign Influenza A/B Antigen Rapid Test is a rapid (result in 15 minutest) lateral flow immunoassay which can be used as an aid for the qualitative detection and differential diagnosis of influenza A and influenza B from anterior nasal swab specimens.

【TEST PRINCIPLE】

The GeneSign Influenza A/B Antigen Rapid Test applies the chemical extraction of viral antigens and the double antibody sandwich method to detect the nucleocapsid protein antigen from influenza A virus and Influenza B virus in anterior nasal swab specimens. The test consists of test lines of monoclonal antibody and pads containing monoclonal antibody-dye conjugate. If the specimen contains sufficient antigen, both the test line (T) and the control line (C) will appear, and the result will be positive. If the specimen does not contain the antigen or no antigen is detected or the level of antigen is below the limit of detection, the test line (T) will not appear, only control line (C) will appear.

When performing the test, an anterior nasal swab specimen is collected and placed into the extraction solution prefilled in extraction tube and manually squeezed, when the antigen is extracted from disrupted virus particles. Take out the test card from sealed foiled pouch and place on a clean, dry and level surface. Add the extracted specimens to the specimen well (S) on the test card, and the solution of extracted specimen will flow to the test strip and migrate through the pads and nitrocellulose membrane of the test strip. The pads contain detector antibodies conjugated to gold dye and the nitrocellulose membrane contains immobilized capture antibodies. If influenza A and/or influenza B antigens are present in the specimen, they will react with anti-influenza A and/or anti-influenza B antibody coupled to gold dye particles, migrate through the membrane as antigen-antibody-dye complexes, bind to the immobilized capture antibody line(s) on the membrane, and generate a colored line in the specific test line region (T), respectively. The rest of the sample and unbound/bound dye complexes continue to migrate to the control line region (C), where immobilized antibodies to the anti-influenza A and/or B antibodies capture the dye complexes and form the control line.

An internal quality control is included in the test, in the form of a colored line appearing in the control line region (C), indicating that the test is functional, and proper and sufficient volume of specimen has been applied to enable migration through the test and control lines, regardless of whether there is a test line or not. If the control line (C) does not appear within the testing time, test result is invalid and the test should be repeated with a new test card and specimen.

【MATERIALS PROVIDED】

- Test card individually foil pouched with a desiccant
- Extraction solution
- Sterile swab
- Instructions for use

【MATERIALS REQUIRED BUT NOT PROVIDED】

Timer, Pipette, Positive Control Swab, Negative Control Swab, Personal Protective Equipment per local recommendations (such as face mask, face shield/eye goggles and gloves), Biohazard Container.

【WARNINGS AND PRECAUTIONS】

1. For *in vitro* diagnostic use only.
2. Do not reuse the test.
3. Do not freeze the test kit or its components.
4. These instructions must be carefully read and strictly followed by a trained healthcare professional to achieve accurate results. All users have to read the instructions before performing test.
5. The test kit is only for the detection of influenza A and B virus antigen, not for any other viruses or pathogens.
6. The test is intended for use with direct anterior nasal swab specimens.
7. Fresh samples are recommended for use to ensure optimal performance. Freshly collected specimens should be tested immediately.
8. Inadequate or inappropriate specimen collection, storage, and transportation are likely to result in false negative test results.
9. Training in specimen collection is highly recommended because of the importance of specimen quality.
10. To collect anterior nasal swab specimens, use the swab provided in this test kit only. Other swabs may not work properly.
11. Do not eat, drink or smoke in the area where handling specimens or performing the test.
12. Do not use the test kit beyond its expiration date.
13. Do not mix components from different kit lots.
14. Leave test card sealed in its foil pouch until just before use. Do not use the test card if the pouch is damaged or the seal is broken.
15. To avoid contamination or inaccurate test result, do not touch the absorbent tip of swab or reaction area of test card when performing the test.
16. Wear appropriate personal protection equipment and gloves when performing the test, collecting and handling patient specimens.
17. Change gloves between handling of suspected specimens and performing each test.
18. Dispose of all used test devices and potentially contaminated materials in a biohazard container as if they were infectious waste and dispose according to applicable local laws and regulations.


【STORAGE AND STABILITY】

1. The test kit should be stored at a temperature between 2-30°C, away from direct sunlight. Do not freeze the kit or its components.
2. The shelf life of the kit is as indicated on the outer package (24 months from date of manufacture).
3. This test kit is stable until the expiration date marked on the outer package and foil pouch. Ensure all test components are at room temperature (15-30°C) before use.
4. Perform the test immediately after taking out the test card from the foil pouch.

【SAMPLE COLLECTION AND STORAGE】

1. Specimen Collection

To achieve accurate test result, good sample collection is the most important first step. Therefore, carefully follow the instructions below to collect anterior nasal swab specimens to obtain as much secretion as possible.

Use the sterile swab provided in the GeneSign Influenza A/B Antigen Rapid Test kit only. Allow the individual's head to relax naturally. Insert the entire soft end of the swab into the individual's nostril no more than 2.5 cm. Slowly rotate the swab in a circular path and gently press against the nostril wall into the nostril of the individual to the nasal palate at least 5 times for no less than 15 seconds. Get as much secretion as possible on the soft end of the swab. Gently withdraw the swab. Repeat the same process in the other nostril using the same swab.

2. Specimen Storage

If the specimen cannot be tested immediately after collection, properly store the specimen in a Viral Transport Media (VTM) or Universal Transport Media (UTM) contained in a lidded storage container. The specimens contained in VTM or UTM can be stored for up to 72 hours when refrigerated (2-8°C) or frozen (-20°C).

【TEST PREPARATION】


Before testing, open the package and equilibrate the test card, extraction solution and specimens to room temperature, and shake the extraction solution gently before use. The most suitable temperature condition to perform the test is room temperature (15-30°C). If the test kit is stored at room temperature, it can be opened and used immediately.

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1. Tear off the sealing film of the extraction solution tube.
2. Place the swab in the extraction solution tube, rotate the swab for about 10 times and no less than 15 seconds, and press the swab tip against the tube wall to release the antigen in the swab.
3. Squeeze the swab over the swab tip to withdraw the swab so as to extract as much liquid as possible from the swab. Dispose of used swabs according to biohazard waste disposal method and local regulations.
4. Install the dripper of the tube and shake well the tube.
5. Take out the test card from sealed foil pouch and place on a dry, clean and level surface. Add two drops of extracted specimen to the specimen well(S) on the test card, and start the timer.

Notes: Applying sufficient amount of extracted specimen is essential for a valid test result. If migration (the wicking of membrane) is not observed in the test window after one minute, add one more drop to specimen well.

6. Read the results at 15 minutes, and the result after 30 minutes is no longer valid.

(Please refer to the illustrations below)

Positive: Presence of both a reddish purple quality control line (C) and a reddish purple test line (T) indicate that influenza A and/or influenza B antigen has been detected. A positive result does not rule out co-infections with other pathogens or identify any specific influenza A virus subtype.


NOTE: The presence of any test line (T), no matter how faint, within the designated observation time, indicates a positive result.

Negative: If there is only a reddish purple quality control line (C) without test line at the T positions, indicates that influenza A and influenza B antigen has not been detected. A negative result does not exclude influenza virus infection.

NOTE: Negative results are presumptive and may need to be confirmed with a molecular assay.

Invalid: If the quality control line (C) is not observed in 15 minutes, the test result will be invalid regardless of whether there is a test line or not, and the test should be repeated using a new test card.

NOTE: Co-infection with Influenza A and B is rare. If test results are positive for both antigens, the specimens should be re-tested.

【LIMITATIONS】

- The test kit is only used for qualitative testing and does not indicate the number of influenza A or B virus antigen in the specimens.
- Failure to follow the instructions or interpretation of test results may adversely affect test performance and/or invalidate the test results.
- A negative test result may occur if the level of antigen in a specimen is below the detection limit of the test.
- Positive test results do not rule out co-infections with other pathogens.
- Positive test results do not identify specific influenza A virus subtypes.
- The test kit can identify influenza A and B viruses, but it cannot differentiate influenza subtypes.
- If differentiation of specific influenza A or B subtypes and strains is needed, it is required to perform additional testing in consultation with state or local public health authorities.
- Individuals who received nasally administered influenza A vaccine may produce positive test results for up to three days after vaccination.
- The performance of this assay has not been evaluated for use in patients without signs and symptoms of respiratory infection.
- The performance of this test has not been evaluated for immunocompromised individuals.
- Negative results do not rule out infection with influenza A, or influenza B and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions.
- Negative results for patients with symptom onset for more than five days should be treated as presumptive and confirmation with a molecular assay, and patient management may be performed according to state or local regulations as needed.
- Positive and negative predictive values are highly dependent on prevalence. False negative test results are more likely during peak activity when prevalence of disease is high. False positive test results are more likely during periods of low activity when prevalence is moderate to low.

【PERFORMANCE CHARACTERISTIC】

1. Analytical sensitivity (Limit of detection, LoD)

To test the limit of detection, heat inactivated viruses of the virus isolates of influenza virus type A and influenza virus type B were used. When tests were repeated 10 times with the diluted virus isolates, the lowest concentration at which more than 95% of the total tested number could be judged as positive test result was determined as the LoD. The LoD of the virus isolates are as follows.

Virus Type	Virus Strains	LoD (TCID ₅₀ /mL)
Influenza A	Liao Ning/1183/2007 (H1N1)	2.2 x 10 ³
Influenza A	Brisbane /10/2007 (H3N2)	3.9 x 10 ³
Influenza A	A/Taiwan/42/06	1.3 x 10 ³
Influenza A	A/HongKong/8/68	2.58 x 10 ⁴
Influenza A	A/Victoria/3/75	1.16 x 10 ²
Influenza A	A/Beijing/302/54	3.6 x 10 ³
Influenza A	A/swine/Guangdong/2/01	6.29 x 10 ³
Influenza A	S-OIV A/HK/415742/09	1.68 x 10 ³
Influenza A	S-OIV A/California/4/09	1.9 x 10 ³
Influenza B	Guang Dong/1512/2010	1.6 x 10 ³
Influenza B	Jiang Xi/32/2000	2.9 x 10 ²
Influenza B	B/ Malaysia/2506/2004	3.2 x 10 ⁵
Influenza B	B/1715	6.7 x 10 ³
Influenza B	B/1704	6.8 x 10 ²
Influenza B	B/Taiwan/2/62	1.58 x 10 ³
Influenza B	B/179	1.98 x 10 ³
Influenza B	B/668	1.5 x 10 ⁴

2. Cross-reactivity

The cross-reactivity against various microorganisms and viruses that may exist in the specimens was tested. The GeneSign Influenza A/B Antigen Rapid Test does not show cross reactivity with the following pathogens up to the concentrations listed in the table below.

Interfering Substances	Concentration	Interfering Substances	Concentration
Staphylococcus aureus	10 ⁶ cfu/mL	Rhinovirus	10 ⁶ pfu/mL
Staphylococcus epidermidis	10 ⁶ cfu/mL	OC43	10 ⁶ pfu/mL
Streptococcus pneumoniae	10 ⁶ cfu/mL	NL63	10 ⁶ pfu/mL
Streptococcus pyogenes	10 ⁶ cfu/mL	229E	10 ⁶ pfu/mL
Candida albicans	10 ⁶ cfu/mL	Adenovirus C1	10 ⁶ pfu/mL
Bordetella pertussis	10 ⁶ cfu/mL	Adenovirus 71	10 ⁶ pfu/mL
Mycoplasma pneumoniae	10 ⁶ cfu/mL	Human Metapneumovirus (hMPV)	10 ⁶ pfu/mL
Chlamydia pneumoniae	10 ⁶ cfu/mL	Parainfluenza virus 1	10 ⁶ pfu/mL
Legionella pneumophila	10 ⁶ cfu/mL	Parainfluenza virus 2	10 ⁶ pfu/mL
Haemophilus influenzae	10 ⁶ cfu/mL	Parainfluenza virus 3	10 ⁶ pfu/mL
Enterovirus	10 ⁶ pfu/mL	Parainfluenza virus 4	10 ⁶ pfu/mL
Respiratory syncytial virus	10 ⁶ pfu/mL	MERS-coronavirus	10 ⁶ pfu/mL
SARS-coronavirus (SARS-CoV)	10 ⁶ pfu/mL	SARS-COV-2	10 ⁶ pfu/mL
Human coronavirus HKU1	10 ⁶ pfu/mL		

3. Endogenous Interference Studies

The potential interference of endogenous substances with the antibodies used for the detection of influenza virus type A and influenza virus type B was examined by testing below substances in a negative clinical matrix, in the absence or presence of each virus, at 3 x LOD concentrations for influenza virus type A and influenza virus type B. The interference study was conducted using medically relevant concentrations of the potentially interfering substances listed below to assess the potential interference of the substances on the performance of the GeneSign Influenza A/B Antigen Rapid Test. No interference was seen with the listed substances when tested at the concentration presented in the table below.

Interfering substances	Concentration	Result
Whole Blood	0.04	No Interference
Mucin	0.005	No Interference
Chloraseptic (Menthol/Benzocaine)	1.5 mg/mL	No Interference
Naso GEL (NeilMed)	5% v/v	No Interference
CVS Nasal Drops (Phenylephrine)	15% v/v	No Interference
Afrin (Oxymetazoline)	15% v/v	No Interference
CVS Nasal Spray (Cromolyn)	15% v/v	No Interference
Zicam	5% v/v	No Interference
Homeopathic (Alkalol)	1:10 dilution	No Interference
Sore Throat Phenol Spray	15% v/v	No Interference
Tobramycin	4 µg/mL	No Interference
Mupirocin	10 mg/mL	No Interference
Fluticasone Propionate	5% v/v	No Interference
Tamiflu (Oseltamivir Phosphate)	5 mg/mL	No Interference

4. High Dose Hook Effect

A high-dose hook effect was not detected in the GeneSign Influenza A/B Antigen Rapid Test for the influenza A and influenza B virus at the concentration listed below.

Virus Type	Virus Strains	Concentration
Influenza A (H1N1)	Liao Ning/1183/2007	1.3 x 10 ⁵ TCID50/mL
Influenza A (H3N2)	Brisbane /10/2007	3.1 x 10 ⁵ TCID50/mL
Influenza B	Guang Dong/1512/2010	1.06 x 10 ⁵ TCID50/mL
Influenza B	Jiang Xi/32/2000	2.8 x 10 ⁵ TCID50/mL

【CLINICAL PERFORMANCE】

To evaluate the clinical performance of the GeneSign Influenza A/B Antigen Rapid Test, directly collected anterior nasal swab specimens which were confirmed with reference RT-PCR assay were tested, including 1181 specimens (163 positive nasal swab specimens, 990 negative nasal swab specimens) for Influenza type A, and 1272 specimens (195 positive nasal swab specimens, 1050 negative nasal swab specimens) for Influenza type B were tested. The positive percent agreement and negative percent agreement are as follows.

Clinical performance compared to reference RT-PCR Assay: Influenza A

GeneSign Influenza A/B Antigen Rapid Test	Reference RT-PCR Assay		
	Positive	Negative	Total
Influenza A Positive	163	16	179
Influenza A Negative	12	990	1002
Total	175	1006	1181

Sensitivity (Positive Percent Agreement): $93.14\% = 163/175$ (95% CI: 88.40%~96.03%)

Specificity (Negative Percent Agreement): $98.40\% = 990/1006$ (95% CI: 97.43%~99.02%)

Accuracy (Overall Percent Agreement): $97.62\% = (163+990)/1181$ (95% CI: 96.59%~98.35%)

Clinical performance compared to reference RT-PCR Assay: Influenza B

GeneSign Influenza A/B Antigen Rapid Test	Reference RT-PCR Assay		
	Positive	Negative	Total
Influenza B Positive	195	15	210
Influenza B Negative	12	1050	1062
Total	207	1065	1272

Sensitivity (Positive Percent Agreement): $94.20\% = 195/207$ (95% CI: 90.14%~96.65%)

Specificity (Negative Percent Agreement): $98.59\% = 1050/1065$ (95% CI: 97.69%~99.14%)

Accuracy (Overall Percent Agreement): $97.87\% = (195+1050)/1272$ (95% CI: 96.93%~98.54%)

【REFERENCES】

- Shaw MW, Arden NH and Massab HF. New aspects of influenza viruses. Clin. Microbiol. Rev. 5: 74-92 (1992)
- WHO recommendations on the use of rapid testing for influenza diagnosis, July 2005.
- Peaper DR, Landry ML. Rapid diagnosis of influenza: state of the art. Clin Lab Med. 2014; 34(2):365-385.doi:10.1016 /j.cll. 2014. 02. 009.

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	Consult instruction for use		For <i>in vitro</i> diagnostic use only		Catalog number		Temperature limit
	Lot number		Use by		Do not reuse		Contains sufficient for <X> tests
	Keep dry		Manufacturer		Date of manufacture		Keep away from sunlight
	Do not use if package is damaged		CE mark		Importer		Distributor
	Authorized Representative in the European Community		Country of manufacture		Model number		Unique device identifier

GeneSign Biotech (Xiamen) Co., Ltd.

Unit 03, 8th Floor, Building B14, Xiamen Biomedical Industrial Park
2076 Wengjiao West Road, Haicang District, Xiamen, Fujian, China
www.genesignbio.com

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SUNGO Europe B.V.