https://trungtamthuoc.com/

2026.04.15 16:
54:02+07'00'

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Email: apac.complaint@alcon.com

Trang web: https://www.alcon.com/contact-us

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2:07 pm, Oct 14, 2025

Alcon PACKAGING

FILE NAME: 300072072-0625 I SYSHYD 10ml VN		
ARTIST: Aranda, Francisco	ART REV: D	DATE: 14.10.2025
SPECIAL INSTRUCTIONS: SPECIAL INSTRUCTIONS GO HERE		

VENDORS PLEASE READ BELOW

COATING = 38 – 60 gloss units as measured by BYK-Gardner Micro-Gloss 60° gloss meter or similar device.
HIGH GLOSS COATING = No less than 75 gloss units as measured by BYK-Gardner Micro-Gloss 60° gloss meter or similar device.
MATTE COATING = No greater than 25 gloss units as measured by BYK-Gardner Micro-Gloss 60° gloss meter or similar device.

INSTRUCTIONS FOR USE

en This Instructions For Use (IFU) leaflet contains important product use and safety information. Please read carefully and keep this for future use

PRODUCT NAME:
SYSTANE™ HYDRATION Lubricant Eye Drops

DESCRIPTION:
SYSTANE™ HYDRATION Lubricant Eye Drops is a sterile solution containing sodium hyaluronate, polyethylene glycol 400, propylene glycol, hydroxypropyl guar, sorbitol, aminomethyl propanol, boric acid, sodium borate, edetate disodium, sodium citrate, potassium chloride, sodium chloride, POLYQUAD™ (polyquaternium-1) 0.001% preservative, and purified water. May contain hydrochloric acid and/or sodium hydroxide to adjust pH.

CLINICAL BENEFITS:
SYSTANE™ HYDRATION Lubricant Eye Drops lubricates the eye surface to provide dry eye symptom relief, protection from further irritation due to dry eye and helps to keep the tear film stable. Additionally, SYSTANE™ HYDRATION Lubricant Eye Drops provides temporary relief from dryness related to contact lens wear.

INTENDED PURPOSE/USE:
SYSTANE™ HYDRATION Lubricant Eye Drops is intended to lubricate the eye surface and rewet soft contact lenses (including silicone hydrogel lenses).

INTENDED USERS/TARGET POPULATION:
Adults (18 years and older) experiencing dry eye symptoms, including contact lens wearers, who are capable of understanding and following the directions and instructions for use (or have a caregiver who can act on their behalf).

INDICATION (Uses):
SYSTANE™ HYDRATION Lubricant Eye Drops is indicated for the temporary relief of burning and irritation in persons experiencing dry eye symptoms.

PERFORMANCE CHARACTERISTICS:
Based on the total body of evidence from clinical investigations, published clinical literature and market experience (post market surveillance), SYSTANE™ HYDRATION Lubricant Eye Drops perform as indicated for the temporary relief of burning and irritation in persons experiencing dry eye symptoms.

DIRECTIONS:

- **Shake well before using.**
- SYSTANE™ HYDRATION Lubricant Eye Drops may be used as needed throughout the day to relieve the symptoms of dry eye.
- Instill 1 or 2 drops in the affected eye(s) as needed and blink.
- While wearing contact lenses, if minor irritation, discomfort, or blurring occurs, instill 1 or 2 drops in the affected eye(s) and blink 2 or 3 times.

CONTRAINDICATIONS (Reasons not to use):
People allergic to any ingredient in SYSTANE™ HYDRATION Lubricant Eye Drops should not use this product.

WARNINGS:
It is essential that you follow your eye care professional's directions and all labeling instructions for proper use of the product. If you experience persistent eye discomfort, eye pain, excessive tearing, vision changes, or redness of the eye, stop using SYSTANE™ HYDRATION Lubricant Eye Drops and consult your eye care professional as the problem could become more serious.

SYSTANE™ HYDRATION Lubricant Eye Drops is not indicated for use as a cleanser or disinfectant for contact lenses. Do not soak or disinfect your lenses in SYSTANE™ HYDRATION Lubricant Eye Drops.

Temporary blurred vision or other visual disturbances may affect the ability to drive or use machines. If blurred vision occurs after application, wait until the vision clears before driving or using machinery.

PRECAUTIONS:

- For ocular use only.
- **TAMPER EVIDENT:** Do not use if security seal is damaged or missing or tamper evident cap is damaged.
- Do not use if package or bottle (sterile packaging) is damaged.

- To avoid contamination, never touch the dropper tip of the container to any surface.
- If swallowed, get medical help.
- If solution changes color or becomes cloudy, do not use.
- Replace cap after using.
- Keep container tightly closed when not in use.
- This product contains a level of boric acid and sodium borate that are safe when used according to label instructions. Boric acid and sodium borate have the potential to cause reproductive toxicity.
- Use before the expiration date marked on the product.
- Once the bottle is opened, discard any remaining solution as instructed below:
For 3 mL fill size, discard remaining solution after three months.
For 5mL and 10 mL fill size, discard remaining solution after six months.
- Store between 2°C to 30°C.
- Keep out of the reach of children.

RESIDUAL RISKS AND UNDESIRABLE SIDE EFFECTS:
Possible residual risks and/or undesirable side effects for SYSTANE™ HYDRATION Lubricant Eye Drops include: allergic response/hypersensitivity, decreased vision (temporary and permanent), eye pain, inflammation, irritation, microbial infection, redness of the eye, systemic toxicity, and toxicity response. If you experience persistent eye discomfort, eye pain, excessive tearing, vision changes or redness of the eye, stop using SYSTANE™ HYDRATION Lubricant Eye Drops and consult your eye care professional as the problem could become more serious.

Serious Incident Reporting:
Any serious incident related to the use of this medical device should be reported to Alcon Laboratories, Inc.:

Contact the local country office or your Alcon distributor.

Email: apac.complaint@Alcon.com
Website: https://www.alcon.com/contact-us

These serious incidents should also be reported to the competent authority for medical devices of your country.

HOW SUPPLIED:
SYSTANE™ HYDRATION Lubricant Eye Drops is supplied in plastic bottles containing fill size of 3 mL or 5 mL or 10 mL.

STERILE A Sterilized using aseptic processing techniques.

GLOSSARY TABLE FOR SYMBOLS	
Symbol	Definitions
	Caution
	Consult instructions for use or consult electronic instructions for use
	Do not use if package is damaged and consult instructions for use
	Date of manufacture
	Contains hazardous substances
	Single sterile barrier system
	Medical Device

WARRANTY:
Not applicable.

STORAGE CONDITION:
Store between 2°C to 30°C.

LEGAL MANUFACTURER:

Alcon Laboratories, Inc.
6201 South Freeway
Fort Worth, Texas 76134, USA

Based on EU Master Labeling Content -
CLEAR ID: TXT000963, Rev 0725

ALC-DFU Systane Hydration-V0825