

Status: Currently Official on 14-Feb-2025  
 Official Date: Official as of 01-Dec-2021  
 Document Type: Reagents  
 DocId: GUID-B7761474-B8FC-4FE9-815B-48F8AEC9845E\_2\_en-US  
 DOI: [https://doi.org/10.31003/USPNF\\_R8085\\_02\\_01](https://doi.org/10.31003/USPNF_R8085_02_01)  
 DOI Ref: 1o2dd

© 2025 USPC  
 Do not distribute

**Add the following:**

## ▲0.01 N Ceric Ammonium Sulfate VS

Transfer 6.5 g of [ceric ammonium sulfate](#) to a 1000-mL volumetric flask containing a mixture of 30 mL of [sulfuric acid](#) and 500 mL of [water](#). Swirl until dissolved. Cool and dilute with [water](#) to volume.

### Standardization

See [Volumetric Solutions, 1. Introduction](#).

See [Titrimetry \(541\)](#).

Standardize by one of the following procedures. [NOTE—Other standardization procedures may be used. See [Volumetric Solutions, 2. Preparation and Standardization, 2.3 Standardization](#).]

**Standardization with visual endpoint:** Transfer 25 mL of the ceric ammonium sulfate solution into an Erlenmeyer flask, add 2 g of [potassium iodide](#) and about 150 mL of [water](#). Mix well. Add 1 mL of [starch TS](#) and titrate immediately with freshly standardized [0.01 M sodium thiosulfate VS](#) until the blue color disappears.

$$N = \frac{\text{mL Na}_2\text{S}_2\text{O}_3 \times N \text{ Na}_2\text{S}_2\text{O}_3}{\text{mL Ce}(\text{NH}_4)_4(\text{SO}_4)_4}$$

**Standardization with potentiometric endpoint:** Transfer 10 mL of the ceric ammonium sulfate solution into an Erlenmeyer flask, add 2 g of [potassium iodide](#) and about 150 mL of [water](#). Mix well. Titrate immediately with freshly standardized [0.01 M sodium thiosulfate VS](#) using a combined platinum electrode.

$$N = \frac{\text{mL Na}_2\text{S}_2\text{O}_3 \times N \text{ Na}_2\text{S}_2\text{O}_3}{\text{mL Ce}(\text{NH}_4)_4(\text{SO}_4)_4}$$

[NOTE—If this volumetric solution is used in a qualitative application such as pH adjustment, dissolution medium, or diluent, its standardization is not required.]

▲ (USP 1-Dec-2021)

**Auxiliary Information** - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
0.01 N CERIC AMMONIUM SULFATE VS	<a href="#">Margareth R.C. Marques</a> Principal Scientific Liaison	HDQ Headquarters

### Most Recently Appeared In:

Pharmacopeial Forum: Volume No. 46(4)

**Current DocID:** [GUID-B7761474-B8FC-4FE9-815B-48F8AEC9845E\\_2\\_en-US](#)

**DOI:** [https://doi.org/10.31003/USPNF\\_R8085\\_02\\_01](https://doi.org/10.31003/USPNF_R8085_02_01)

**DOI ref:** [1o2dd](#)