

Status: Currently Official on 14-Feb-2025
 Official Date: Official as of 01-May-2021
 Document Type: Reagents
 DocId: GUID-D0F17316-A5C8-43DD-994E-5150A9998414_2_en-US
 DOI: https://doi.org/10.31003/USPNF_R3370_02_01
 DOI Ref: mgh0w

© 2025 USPC
 Do not distribute

Change to read:

0.05 M Edetate Disodium VS

▲ (USP 1-May-2021)

Dissolve 18.6 g of [edetate disodium](#) in [water](#) to make 1000 mL.

Standardization

▲ See [Volumetric Solutions, 1. Introduction](#).

See [Titrimetry \(541\)](#).

Standardize by one of the following procedures. [NOTE—Other standardization procedures may be used. See [Volumetric Solutions, 2. Preparation and Standardization, 2.3 Standardization](#).]

Standardization with visual endpoint: ▲ (USP 1-May-2021) Accurately weigh about 200 mg of [chelometric standard calcium carbonate](#),

previously dried according to the label instructions or, if this information is not available, at ▲210° for 4 h, ▲ (USP 1-May-2021) transfer to a 400-mL beaker, add 10 mL of [water](#), and swirl to form a slurry. Cover the beaker with a watch glass, and introduce 2 mL of [diluted hydrochloric acid](#) from a pipet inserted between the lip of the beaker and the edge of the watch glass. Swirl the contents of the beaker to dissolve the calcium carbonate. Wash down the sides of the beaker, the outer surface of the pipet, and the watch glass with [water](#), and dilute with [water](#) to about 100 mL. While stirring the solution, preferably with a magnetic stirrer, add about 30 mL of the edetate disodium solution from a 50-mL buret. ▲ Adjust the solution with [sodium hydroxide TS](#) to a pH of 12–13, add ▲ (USP 1-May-2021) 300 mg of [hydroxy naphthol blue](#), and continue the titration with the edetate disodium solution to a blue endpoint. ▲

$$M = \frac{(gCaCO_3) \times (Assay / 100) \times 1000}{100.09 \times mL EDTA}$$

Standardization with potentiometric endpoint: Accurately weigh about 50 mg of [chelometric standard calcium carbonate](#), previously dried according to the label instructions or, if this information is not available, at 210° for 4 h, transfer to a 400-mL beaker, add 10 mL of [water](#), and swirl to form a slurry. Cover the beaker with a watch glass, and introduce 0.5 mL of [diluted hydrochloric acid](#) from a pipet inserted between the lip of the beaker and the edge of the watch glass. Swirl the contents of the beaker to dissolve the calcium carbonate. After dissolution, wash the sides of the beaker, the outer surface of the pipet, and the watch glass with [water](#), and dilute with [water](#) to about 100 mL. Add 4 mL of [sodium hydroxide TS](#) and titrate with the edetate disodium solution using a combined calcium ion-selective electrode.

$$M = \frac{(gCaCO_3) \times (Assay / 100) \times 1000}{100.09 \times mL EDTA}$$

▲ (USP 1-May-2021)

[NOTE—If this volumetric solution is used in a qualitative application such as pH adjustment, dissolution medium, or diluent, its standardization is not required.]

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
0.05 M EDETATE DISODIUM VS	Margareth R.C. Marques Principal Scientific Liaison	HDQ Headquarters

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. 45(6)

Current DocID: [GUID-D0F17316-A5C8-43DD-994E-5150A9998414_2_en-US](#)

DOI: https://doi.org/10.31003/USPNF_R3370_02_01

DOI ref: [mgh0w](#)