

Status: Currently Official on 13-Feb-2025
Official Date: Official Prior to 2013
Document Type: Reagents
DocId: GUID-13524F14-09D3-4DF4-81A9-F576A89E9C0C_1_en-US
DOI: https://doi.org/10.31003/USPNF_R3069_01_01
DOI Ref: jr9gz

© 2025 USPC
Do not distribute

Acetate Buffer TS

—Dissolve 320 g of ammonium acetate in 500 mL of water, add 5 mL of glacial acetic acid, dilute with water to 1000.0 mL, and mix. This solution has a pH between 5.9 and 6.0.

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
ACETATE BUFFER TS	Margareth R.C. Marques Principal Scientific Liaison	HDQ Headquarters

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. Information currently unavailable

Current DocID: [GUID-13524F14-09D3-4DF4-81A9-F576A89E9C0C_1_en-US](#)

DOI: https://doi.org/10.31003/USPNF_R3069_01_01

DOI ref: [jr9gz](#)

OFFICIAL