

Status: Currently Official on 13-Feb-2025
Official Date: Official Prior to 2013
Document Type: USP Monographs
DocId: GUID-474A5736-9D2F-4746-A89A-DA2D4E329626_1_en-US
DOI: https://doi.org/10.31003/USPNF_M548_01_01
DOI Ref: 9iv9m

© 2025 USPC
Do not distribute

Acetohydroxamic Acid Tablets

DEFINITION

Acetohydroxamic Acid Tablets contain NLT 90.0% and NMT 110.0% of the labeled amount of acetohydroxamic acid ($C_2H_5NO_2$).

IDENTIFICATION

- **A.** Tablets produce a purple color when mixed with an acidic solution of ferric chloride.

ASSAY

• PROCEDURE

Ferric chloride solution: 20 mg/mL of ferric chloride in 0.1 N hydrochloric acid

Standard solution: 500 µg/mL of [USP Acetohydroxamic Acid RS](#) in 0.1 N hydrochloric acid

Sample solution: Weigh, and finely powder NLT 20 Tablets. Transfer an accurately weighed portion of the powder, equivalent to about 500 mg of acetohydroxamic acid, to a 1000-mL volumetric flask, add about 500 mL of 0.1 N hydrochloric acid, and shake for 1 min. Dilute with 0.1 N hydrochloric acid to volume, and mix. Filter, discarding the first 40 mL of the filtrate. Use the clear filtrate.

Blank: 0.1 N hydrochloric acid

Analysis

Samples: *Standard solutions, Sample solution, and Blank*

Transfer 10.0 mL each of the *Standard solution, Sample solution, and Blank* to separate 100-mL volumetric flasks. To each flask add 50 mL of 0.1 N hydrochloric acid and 10.0 mL of *Ferric chloride solution*, and dilute with 0.1 N hydrochloric acid to volume. Without delay, concomitantly determine the absorbances of the solutions at the wavelength of maximum absorbance at about 502 nm, using the *Blank* to set the instrument.

Calculate the percentage of labeled amount of acetohydroxamic acid ($C_2H_5NO_2$) in the portion of Tablets taken:

$$\text{Result} = (A_U/A_S) \times (C_S/C_U) \times 100$$

A_U = absorbance of the *Sample solution*

A_S = absorbance of the *Standard solution*

C_S = concentration of [USP Acetohydroxamic Acid RS](#) in the *Standard solution* (µg/mL)

C_U = nominal concentration of acetohydroxamic acid in the *Sample solution* (µg/mL)

Acceptance criteria: 90.0%–110.0%

PERFORMANCE TESTS

- [DISSOLUTION, Procedure for a Pooled Sample\(711\).](#)

Medium: 0.01 N hydrochloric acid; 900 mL

Apparatus 1: 100 rpm

Time: 30 min

Analysis: Calculate the percentage of the labeled amount of acetohydroxamic acid ($C_2H_5NO_2$) dissolved, using the procedure in the Assay, using a filtered portion of the solution under test as *Sample solution* in comparison with a *Standard solution* having a known concentration of [USP Acetohydroxamic Acid RS](#) in *Medium*.

Tolerances: NLT 85% (Q) of the labeled amount of acetohydroxamic acid ($C_2H_5NO_2$) is dissolved.

- [UNIFORMITY OF DOSAGE UNITS \(905\):](#) Meet the requirements

IMPURITIES

• LIMIT OF HYDROXYLAMINE

Buffer: 1.36 g/L of monobasic potassium phosphate in water, adjusted with 1 M potassium hydroxide to a pH of 7.4

Solution A: 1 mg/mL of pyridoxal 5-phosphate monohydrate in *Buffer*, prepared in a low-actinic flask fresh before use

Standard stock solution: 2.0 mg/mL of hydroxylamine hydrochloride in water

Standard solutions: Transfer 5.0, 10.0, and 15.0 mL of the *Standard stock solution* to separate 100-mL volumetric flasks, and dilute with water to volume.

Sample solution: Weigh, and finely powder NLT 20 Tablets. Transfer a portion of the powder, equivalent to about 1500 mg of acetohydroxamic acid to a 50-mL stoppered centrifuge tube. Add 30.0 mL of water, shake for about 2 min, and centrifuge. Pipet 15.0 mL of the clear solution into a 50-mL beaker, add just enough water to cover the electrode of a calibrated pH meter, and while stirring, adjust with 0.5 M potassium hydroxide to a pH of 7.4. Quantitatively transfer the contents of the beaker with the aid of small portions of water to a 50-mL volumetric flask, dilute with water to volume, and mix.

Blank: Water

Analysis

Samples: *Standard solutions, Sample solution, and Blank*

Transfer 2.0 mL of each *Standard solution*, the *Sample solution*, and *Blank* into separate 100-mL volumetric flasks. To each flask add 4.0 mL of *Solution A*. After 8 min, accurately timed, dilute the contents of each flask with *Buffer* to volume.

Immediately determine the fluorescence intensities of the solutions from the *Standard solutions* and the *Sample solution* in a fluorometer at an excitation wavelength of 350 nm and an emission wavelength of 450 nm, setting the instrument to zero with the *Blank*. Determine the best-fit straight line from the fluorescence intensities of the three *Standard solutions* versus the hydroxylamine hydrochloride concentrations, in µg/mL. From the best-fit straight line, determine the concentration, in µg/mL, of hydroxylamine hydrochloride in the *Sample solution*.

Calculate the percentage of hydroxylamine in the portion of Tablets taken:

$$\text{Result} = (C_U/C) \times (M_{r1}/M_{r2}) \times 100$$

C_U = concentration of hydroxylamine hydrochloride in the *Sample solution* (µg/mL)

C = nominal concentration of acetohydroxamic acid in the *Sample solution* (µg/mL)

M_{r1} = molecular weight of hydroxylamine, 33.03

M_{r2} = molecular weight of hydroxylamine hydrochloride, 69.50

Acceptance criteria: NMT 0.5%

ADDITIONAL REQUIREMENTS

- **PACKAGING AND STORAGE:** Preserve in tight containers.
- **USP REFERENCE STANDARDS** (11).
[USP Acetohydroxamic Acid RS](#)

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
ACETOHYDROXAMIC ACID TABLETS	Documentary Standards Support	SM32020 Small Molecules 3

Chromatographic Database Information: [Chromatographic Database](#)

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. PF 30(1)

Current DocID: GUID-474A5736-9D2F-4746-A89A-DA2D4E329626_1_en-US

DOI: https://doi.org/10.31003/USPNF_M548_01_01

DOI ref: [9iv9m](#)