Status: Currently Official on 13-Feb-2025
Official Date: Official as of 01-Jan-2018
Document Type: USP Monographs
DocId: GUID-276C2A73-5372-482B-B9BD-005237C5EB9A_3_en-US
DOI: https://doi.org/10.31003/USPNF_M2040_03_01
DOI Ref: 8eg09

© 2025 USPC Do not distribute

Aluminum Chloride

 $AICI_{3} \cdot 6H_{2}O$ 241.43 $AICI_{2}$ 133.34

Aluminum chloride, hexahydrate;

Aluminum chloride hexahydrate CAS RN®: 7784-13-6; UNII: 3CYT62D3GA.

Anhydrous CAS RN®: 7446-70-0; UNII: LIF1N9568Y.

DEFINITION

Aluminum Chloride contains NLT 95.0% and NMT 102.0% of aluminum chloride (AlCl₂), calculated on the anhydrous basis.

IDENTIFICATION

• A. IDENTIFICATION TESTS—GENERAL, Aluminum(191) and Chloride(191)

Sample solution: 100 mg/mL

Acceptance criteria: Meets the requirements

ASSAY

• PROCEDURE

Edetate disodium titrant: Prepare and standardize as directed in Reagents, Volumetric Solutions, Edetate Disodium, Twentieth-Molar (0.05 M).

Sample solution: 20 mg/mL of aluminum chloride in water

Titrimetric system

Mode: Back-titration

Titrant: 0.05 M zinc sulfate VS **Endpoint detection:** Visual

Analysis: Transfer 10.0 mL of the Sample solution into a 250-mL beaker, and add, in the order named and with continuous stirring, 25.0 mL of Edetate disodium titrant and 20 mL of acetic acid-ammonium acetate buffer TS, and boil gently for 5 min. Cool, and add 50 mL of alcohol and 2 mL of dithizone TS. Titrate the excess edetate disodium with Titrant to a bright rose-pink color. Perform a blank determination, substituting 10 mL of water for the Sample solution, and make any necessary correction. Each mL of Edetate disodium titrant consumed is equivalent to 6.667 mg of aluminum chloride (AICI₂).

Acceptance criteria: 95.0%-102.0% on the anhydrous basis

IMPURITIES

• LIMIT OF SULFATE

Sample solution: 10 mg/mL

Analysis: Add 0.2 mL of barium chloride TS to 10 mL of the Sample solution.

Acceptance criteria: No turbidity is produced within 1 min.

• <u>Iron (241)</u> **Sample:** 1.0 g

Analysis: Dissolve the Sample in 45 mL of water, and add 2 mL of hydrochloric acid.

Acceptance criteria: NMT 10 ppm
• LIMIT OF ALKALIES AND ALKALINE EARTHS

Sample: 1.0 q

Analysis: To a boiling solution of the *Sample* in 150 mL of water add a few drops of methyl red TS, then add 6 N ammonium hydroxide until the color of the solution just changes to a distinct yellow. Add hot water to restore the volume to 150 mL, and filter while hot. Evaporate 75 mL of the filtrate to dryness, and ignite to a constant weight.

Acceptance criteria: 0.5%; the weight of the residue is NMT 2.5 mg.

SPECIFIC TESTS

• Water Determination, Method I (921): 42.0%-48.0%

ADDITIONAL REQUIREMENTS

• Packaging and Storage: Preserve in tight containers.

Auxiliary Information - Please check for your question in the FAQs before contacting USP.

Topic/Question	Contact	Expert Committee
ALUMINUM CHLORIDE	Documentary Standards Support	SM32020 Small Molecules 3

Chromatographic Database Information: Chromatographic Database

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. Information currently unavailable

Current DocID: GUID-276C2A73-5372-482B-B9BD-005237C5EB9A_3_en-US Previous DocID: GUID-276C2A73-5372-482B-B9BD-005237C5EB9A_1_en-US

DOI: https://doi.org/10.31003/USPNF_M2040_03_01

DOI ref: 8eg09