
Status: Currently Official on 17-Feb-2025
Official Date: Official Prior to 2013
Document Type: USP Monographs
Docld: GUID-37049FC7-B8C9-4E71-81F3-BA829B2F8675\_2\_en-US
DOI: https://doi.org/10.31003/USPNF\_M6390\_02\_01
DOI Ref: q2khy

© 2025 USPC Do not distribute

# **Atropine**



C<sub>17</sub>H<sub>23</sub>NO<sub>3</sub> 289.37

Benzeneacetic acid,  $\alpha$ -(hydroxymethyl)-8-methyl-8-azabicyclo[3.2.1]oct-3-yl ester,  $endo-(\pm)$ -;

 $1\alpha H$ ,  $5\alpha H$ -Tropan- $3\alpha$ -ol (±)-tropate (ester) CAS RN<sup>®</sup>: 51-55-8; UNII: 7C0697DR9I.

### **DEFINITION**

Atropine contains NLT 99.0% and NMT 100.5% of atropine (C<sub>17</sub>H<sub>23</sub>NO<sub>3</sub>), calculated on the anhydrous basis.

[CAUTION—Handle Atropine with exceptional care, because it is highly potent.]

## **IDENTIFICATION**

٠Α.

Standard: 36 mg of USP Atropine Sulfate RS

Sample: 30 mg

**Analysis:** Dissolve the *Standard* and *Sample* in individual 60-mL separators with the aid of 5-mL portions of water. To each separator add 1.5 mL of 1 N sodium hydroxide solution and 10 mL of chloroform. Shake for 1 min, allow the layers to separate, and pass the chloroform extracts through separate filters of 2 g of anhydrous granular sodium sulfate supported on pledgets of glass wool. Extract each aqueous layer with two additional 10-mL portions of chloroform, filtering and combining with the respective main extracts. Evaporate the chloroform solutions under reduced pressure to dryness, and dissolve each residue in 10 mL of carbon disulfide.

**Acceptance criteria:** The IR absorption spectrum, determined in a 1-mm cell, of the solution of the *Sample* exhibits maxima only at the same wavelengths as that of the solution of the *Standard*.

• В.

Sample solution: A solution (1 in 50) in 3 N hydrochloric acid

Analysis: Add gold chloride TS to the Sample solution.

Acceptance criteria: A lusterless precipitate is formed (distinction from hyoscyamine, which, similarly treated, yields a lustrous precipitate).

## **ASSAY**

• Procedure

Sample: 400 mg of Atropine

**Analysis:** Dissolve the *Sample* in 50 mL of glacial acetic acid. Titrate with 0.1 N perchloric acid VS to a green endpoint, using 1 drop of crystal violet TS. Perform a blank determination (see <u>Titrimetry (541)</u>). Each mL of 0.1 N perchloric acid is equivalent to 28.94 mg of atropine  $(C_{17}H_{23}NO_2)$ .

Acceptance criteria: 99.0%-100.5% on the anhydrous basis

## **IMPURITIES**

• Residue on Ignition (281): NMT 0.1%

• LIMIT OF FOREIGN ALKALOIDS AND OTHER IMPURITIES

Standard solution: 24 mg/mL of USP Atropine Sulfate RS in methanol

**Sample solution A:** 20 mg/mL of Atropine in methanol **Sample solution B:** 1 mg/mL of Atropine in methanol

**Chromatographic system** 

(See Chromatography (621), Thin-Layer Chromatography.)

Mode: TLC

Adsorbent: 0.5-mm layer of chromatographic silica gel

Application volume: See Analysis.

Developing solvent system: Chloroform, acetone, and diethylamine (5:4:1)

Spray reagent: Potassium iodoplatinate TS

**Analysis** 

**Samples:** Standard solution, 5  $\mu$ L; Sample solution A, 25  $\mu$ L; Sample solution B, 1  $\mu$ L

Apply the Samples to the TLC plate. Allow the spots to dry, and develop the chromatogram in the Developing solvent system until the solvent front has moved three-fourths of the length of the plate. Allow the solvent to evaporate. Locate the spots on the plate by spraying with Spray reagent.

Acceptance criteria: NMT 0.2%; the  $R_F$  value of the principal spot of each Sample solution corresponds to that of the Standard solution; no secondary spot of Sample solution A exhibits intensity equal to or greater than the principal spot of Sample solution B.

• READILY CARBONIZABLE SUBSTANCES TEST (271)

Sample solution: 200 mg in 5 mL of 2 N sulfuric acid

**Acceptance criteria:** The solution has no more color than *Matching Fluid A*, and the solution is colored no more than light yellow upon the addition of 0.2 mL of nitric acid.

#### **SPECIFIC TESTS**

• OPTICAL ROTATION, Angular Rotation (781)

**Sample solution:** 1 g, previously dried at 105° for 1 h, in sufficient 50% alcohol (w/w) to obtain a volume of 20 mL at 25° (using a 200-mm tube)

Acceptance criteria: -0.70° to +0.05° (limit of hyoscyamine)

- Melting Range or Temperature (741): 114°-118°
- Water Determination, Method I(921): NMT 0.2%

### **ADDITIONAL REQUIREMENTS**

- Packaging and Storage: Preserve in tight, light-resistant containers.
- USP REFERENCE STANDARDS (11)

USP Atropine Sulfate RS

Auxiliary Information - Please check for your question in the FAQs before contacting USP.

Topic/Question	Contact	Expert Committee
ATROPINE	Documentary Standards Support	SM42020 Small Molecules 4
REFERENCE STANDARD SUPPORT	RS Technical Services RSTECH@usp.org	SM42020 Small Molecules 4

Chromatographic Database Information: Chromatographic Database

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. PF 43(1)

Current DocID: GUID-37049FC7-B8C9-4E71-81F3-BA829B2F8675\_2\_en-US Previous DocID: GUID-37049FC7-B8C9-4E71-81F3-BA829B2F8675\_1\_en-US

DOI: https://doi.org/10.31003/USPNF\_M6390\_02\_01

DOI ref: q2khy