

Status: Currently Official on 14-Feb-2025
Official Date: Official Prior to 2013
Document Type: USP Monographs
DocId: GUID-E13D12B4-93D3-4087-B463-99CFF4FBD7F1_2_en-US
DOI: https://doi.org/10.31003/USPNF_M14118_02_01
DOI Ref: 3r1vy

© 2025 USPC
Do not distribute

Cefprozil Tablets

DEFINITION

Cefprozil Tablets contain NLT 90.0% and NMT 120.0% of the labeled amount of cefprozil ($C_{18}H_{19}N_3O_5S$).

IDENTIFICATION

• A. THIN-LAYER CHROMATOGRAPHY

Diluent: Acetone and 0.1 N hydrochloric acid (4:1)

Standard solution: 5 mg/mL of [USP Cefprozil \(Z\)-Isomer RS](#) in *Diluent*

Sample solution: Nominally 2.5 mg/mL of cefprozil from 1 Tablet in *Diluent*. Shake for 5 min, and allow the mixture to settle. Use the supernatant.

Chromatographic system

(See [Chromatography \(621\)](#), [Thin-Layer Chromatography](#).)

Adsorbent: 0.25-mm layer of chromatographic silica gel mixture

Application volume: 10 μ L

Developing solvent system: Butyl alcohol, glacial acetic acid, and water (60:20:20)

Analysis

Samples: *Standard solution* and *Sample solution*

Allow the spots to dry, and develop the chromatogram in an equilibrated chamber with the *Developing solvent system* until the solvent front has moved three-fourths of the length of the plate. Remove the plate, and allow it to air-dry in a hood. Place the dry plate in a chamber containing iodine vapors. Examine the plate, and locate the spots.

Acceptance criteria: The R_F value of the principal spot of the *Sample solution* corresponds to that of the *Standard solution*.

• B. The retention times of the cefprozil (Z)-isomer and cefprozil (E)-isomer peaks of the *Sample solution* correspond to those of *Standard solution 1* and *Standard solution 2*, respectively, as obtained in the Assay.

ASSAY

• PROCEDURE

Buffer: 11.5 g/L of monobasic ammonium phosphate in water. Adjust, if necessary, with phosphoric acid to a pH of 4.4.

Mobile phase: Acetonitrile and *Buffer* (100:900). [NOTE—Decreasing the proportion of acetonitrile increases retention times and improves the resolution between the cefprozil isomer peaks.]

Standard solution 1: 0.25 mg/mL of [USP Cefprozil \(Z\)-Isomer RS](#). Use this solution within 6 h.

Standard stock solution: 0.25 mg/mL of [USP Cefprozil \(E\)-Isomer RS](#)

Standard solution 2: 0.025 mg/mL of [USP Cefprozil \(E\)-Isomer RS](#) in water from the *Standard stock solution*. Use this solution within 6 h.

System suitability solution: A mixture of equal volumes of *Standard solution 1* and the *Standard stock solution*. Use this solution within 6 h.

Sample stock solution: Nominally 6 mg/mL of cefprozil in water from Tablets, prepared as follows. Transfer a suitable number of Tablets to a volumetric flask containing water. Allow the Tablets to disintegrate with the aid of swirling and sonication. Dilute with water to volume.

Sample solution: Nominally 0.3 mg/mL of cefprozil in water from the *Sample stock solution*. Pass a portion of this solution through a filter of 0.5- μ m or finer pore size. Use this solution within 6 h.

Chromatographic system

(See [Chromatography \(621\)](#), [System Suitability](#).)

Mode: LC

Detector: UV 280 nm

Column: 3.9-mm \times 25-cm; 5- μ m packing L1

Flow rate: 1 mL/min

Injection volume: 10 μ L

System suitability

Samples: *Standard solution 1* and *System suitability solution*

[NOTE—The relative retention times for cefprozil (Z)-isomer and cefprozil (E)-isomer are about 0.7 and 1.0, respectively.]

Suitability requirements

Resolution: NLT 2.5 between cefprozil (Z)-isomer and cefprozil (E)-isomer, *System suitability solution*

Column efficiency: NLT 2500 theoretical plates, cefprozil (Z)-isomer, *Standard solution 1*

Calculate as follows:

$$\text{Result} = (t_R/W_{h/2})^2 \times 5.545$$

t_R = retention time of cefprozil (Z)-isomer

$W_{h/2}$ = peak width at half-height

Tailing factor: 0.9–1.1, cefprozil (Z)-isomer, *Standard solution 1*

Calculate as follows:

$$\text{Result} = W_{0.1}/2f$$

$W_{0.1}$ = width of the peak at 10% height

f = distance from the peak maximum to the leading edge of the peak measured at 10% of the peak height

Relative standard deviation: NMT 2.0%, *Standard solution 1*

Analysis

Samples: *Standard solution 1*, *Standard solution 2*, and *Sample solution*

Calculate the concentration, in mg/mL, of the cefprozil (Z)-isomer in the *Sample solution*:

$$\text{Result} = (r_U/r_S) \times C_S \times P \times F$$

r_U = peak response of the cefprozil (Z)-isomer from the *Sample solution*

r_S = peak response of the cefprozil (Z)-isomer from *Standard solution 1*

C_S = concentration of [USP Cefprozil \(Z\)-Isomer RS](#) in *Standard solution 1* (mg/mL)

P = potency of the cefprozil (Z)-isomer in [USP Cefprozil \(Z\)-Isomer RS](#) (µg/mg)

F = correction factor, 0.001 mg/µg

Calculate the concentration, in mg/mL, of the cefprozil (E)-isomer in the *Sample solution*:

$$\text{Result} = (r_U/r_S) \times C_S \times P \times F$$

r_U = peak response of the cefprozil (E)-isomer from the *Sample solution*

r_S = peak response of the cefprozil (E)-isomer from *Standard solution 2*

C_S = concentration of [USP Cefprozil \(E\)-Isomer RS](#) in *Standard solution 2* (mg/mL)

P = potency of the cefprozil (E)-isomer in [USP Cefprozil \(E\)-Isomer RS](#) (µg/mg)

F = correction factor, 0.001 mg/µg

Calculate the percentage of the labeled amount of cefprozil ($C_{18}H_{19}N_3O_5S$) in the portion of Tablets taken:

$$\text{Result} = [(C_Z + C_E)/C_U] \times 100$$

C_Z = concentration of the cefprozil (Z)-isomer in the *Sample solution* (mg/mL)

C_E = concentration of the cefprozil (E)-isomer in the *Sample solution* (mg/mL)

C_U = nominal concentration of cefprozil in the *Sample solution* (mg/mL)

Acceptance criteria: 90.0%–120.0%

PERFORMANCE TESTS

• [DISSOLUTION \(711\)](#)

Medium: Water; 900 mL

Apparatus 1: 100 rpm

Time: 45 min

Mobile phase, Standard solution 1, Standard solution 2, System suitability solution, Chromatographic system, and System

suitability: Proceed as directed in the Assay.

Sample solution: Pass the solution under test through a suitable filter. Dilute, if necessary, to 0.3 mg/mL of cefprozil.

Analysis

Samples: *Standard solution 1*, *Standard solution 2*, and *Sample solution*

Calculate the quantity, in mg, of cefprozil (Z)-isomer dissolved:

$$\text{Result} = (r_U/r_S) \times C_S \times D \times V \times P \times F$$

r_U = peak response of cefprozil (Z)-isomer from the *Sample solution*

r_S = peak response of cefprozil (Z)-isomer from *Standard solution 1*

C_S = concentration of [USP Cefprozil \(Z\)-Isomer RS](#) in *Standard solution 1* (mg/mL)

D = dilution factor of the *Sample solution*

V = volume of *Medium*, 900 mL

P = potency of the cefprozil (Z)-isomer in [USP Cefprozil \(Z\)-Isomer RS](#) (µg/mg)

F = correction factor, 0.001 mg/µg

Calculate the quantity, in mg, of cefprozil (E)-isomer dissolved:

$$\text{Result} = (r_U/r_S) \times C_S \times D \times V \times P \times F$$

r_U = peak response of cefprozil (E)-isomer from the *Sample solution*

r_S = peak response of cefprozil (E)-isomer from *Standard solution 2*

C_S = concentration of [USP Cefprozil \(E\)-Isomer RS](#) in *Standard solution 2* (mg/mL)

D = dilution factor of the *Sample solution*

V = volume of *Medium*, 900 mL

P = potency of the cefprozil (E)-isomer in [USP Cefprozil \(E\)-Isomer RS](#) (µg/mg)

F = correction factor, 0.001 mg/µg

Calculate the percentage of cefprozil ($C_{18}H_{19}N_3O_5S$) dissolved:

$$\text{Result} = (M_Z + M_E) \times 100/L$$

M_Z = quantity of cefprozil (Z)-isomer dissolved (mg)

M_E = quantity of cefprozil (E)-isomer dissolved (mg)

L = label claim of cefprozil (mg/Tablet)

Tolerances: NLT 75% (Q) of the labeled amount of cefprozil ($C_{18}H_{19}N_3O_5S$) is dissolved.

- [UNIFORMITY OF DOSAGE UNITS \(905\)](#): Meet the requirements

SPECIFIC TESTS

- [WATER DETERMINATION, Method I\(921\)](#): NMT 7.0%

ADDITIONAL REQUIREMENTS

- **PACKAGING AND STORAGE:** Preserve in tight containers.

- [USP REFERENCE STANDARDS \(11\)](#).

[USP Cefprozil \(E\)-Isomer RS](#)

[USP Cefprozil \(Z\)-Isomer RS](#)

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
CEFPROZIL TABLETS	Documentary Standards Support	SM12020 Small Molecules 1

Chromatographic Database Information: [Chromatographic Database](#)

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. Information currently unavailable

Current DocID: GUID-E13D12B4-93D3-4087-B463-99CFF4FBD7F1_2_en-US

Previous DocID: GUID-E13D12B4-93D3-4087-B463-99CFF4FBD7F1_1_en-US

DOI: https://doi.org/10.31003/USPNF_M14118_02_01

OFFICIAL