Status: Currently Official on 14-Feb-2025
Official Date: Official as of 01-May-2020
Document Type: NF Monographs
Docld: GUID-478339C6-47ED-41ED-AA0C-D0B5316E97FC_2_en-US
DOI: https://doi.org/10.31003/USPNF_M32170_02_01
DOI Ref: n2a5x

© 2025 USPC Do not distribute

Ethylparaben

 $C_9 H_{10} O_3$

166.17

Benzoic acid, 4-hydroxy-, ethyl ester;

Ethyl p-hydroxybenzoate CAS RN[®]: 120-47-8.

DEFINITION

Ethylparaben contains NLT 98.0% and NMT 102.0% of C_oH₁₀O₃.

IDENTIFICATION

Change to read:

- A. Spectroscopic Identification Tests (197), Infrared Spectroscopy: 197M (CN 1-May-2020)
- B. Melting Range or Temperature (741): 115°-118°

ASSAY

• PROCEDURE

Mobile phase, Sample solution, Standard solution B, and **Chromatographic system:** Proceed as described in the procedure for *Related Substances*.

System suitability

Sample: Standard solution B Suitability requirements

Relative standard deviation: NMT 0.85% for 6 injections

Analysis

Samples: Sample solution and Standard solution B

Calculate the percentage of Ethylparaben in the Sample solution:

Result =
$$P \times (r_{II} \times C_{S})/(r_{S} \times C_{II})$$

P = labeled purity of <u>USP Ethylparaben RS</u> expressed as a percentage

r, = peak area of ethylparaben from the Sample solution

C_s = concentration of ethylparaben in Standard solution B

 r_S = peak area of ethylparaben from Standard solution B

C₁₁ = concentration of Ethylparaben in the Sample solution

Acceptance criteria: 98.0%-102.0%

IMPURITIES

Inorganic Impurities

• Residue on Ignition (281): NMT 0.1%, determined on 1.0 g

Organic Impurities

• PROCEDURE: RELATED SUBSTANCES

Mobile phase: Methanol and a 6.8-g/L solution of potassium dihydrogen phosphate (65:35 v/v)

Sample solution: Dissolve 50.0 mg of Ethylparaben in 2.5 mL of methanol, and dilute with *Mobile phase* to 50.0 mL. Dilute 10.0 mL of this solution with *Mobile phase* to 100.0 mL.

Standard solution A: 5.0 µg/mL each of p-hydroxy benzoic acid, <u>USP Methylparaben RS</u>, and <u>USP Ethylparaben RS</u> in Mobile phase

Standard solution B: Dissolve 50.0 mg of <u>USP Ethylparaben RS</u> in 2.5 mL of methanol, and dilute with *Mobile phase* to 50.0 mL. Dilute 10.0 mL of this solution with *Mobile phase* to 100.0 mL.

Standard solution C: Dilute 1.0 mL of the Sample solution with Mobile phase to 20.0 mL. Dilute 1.0 mL of this solution with Mobile phase to 10.0 mL.

Chromatographic system

(See Chromatography (621), System Suitability.)

Mode: LC

Detector: UV 272 nm

Column: 4.6-mm × 15-cm; 5-µm packing L1

Flow rate: 1.3 mL/min Injection size: 10 µL

Run time: About 4 times the retention time of ethylparaben

System suitability

Sample: Standard solution A

[Note—The retention time of ethylparaben is about 3.0 min; the relative retention times for *p*-hydroxybenzoic acid, ethylparaben, and methylparaben are about 0.5, 1.0, and 0.8, respectively.]

Suitability requirements

Resolution: NLT 2.0 between the methylparaben and ethylparaben peaks

Analysis

Samples: Sample solution and Standard solution C

[Note-Disregard any limit that is 0.2 times the area of the principal peak in the chromatogram obtained with Standard solution C (0.1%).]

Acceptance criteria

p-Hydroxybenzoic acid: The peak area in the *Sample solution*, multiplied by 1.4 to correct for the calculation of content, is NMT the area of the principal peak in *Standard solution C* (0.5%).

Unspecified impurities: The peak area of each impurity in the *Sample solution* is NMT the area of the principal peak in *Standard solution C* (0.5%).

Total impurities: The total peak area for all impurities in the *Sample solution* is NMT twice the area of the principal peak in *Standard solution C* (1.0%).

SPECIFIC TESTS

Color of Solution

Sample solution: 100 mg/mL in alcohol

Comparison solution: Mix 2.4 mL of ferric chloride CS, 1.0 mL of cobaltous chloride CS, and 0.4 mL of cupric sulfate CS with 0.3 N hydrochloric acid to make 10 mL. Dilute 5 mL of this solution with 0.3 N hydrochloric acid to make 100 mL. [Note—Prepare and use this solution immediately.]

Analysis

Samples: Alcohol, Sample solution, and Comparison solution

Make the comparison by viewing the solutions downward in matched color-comparison tubes against a white surface (see <u>Color and Achromicity (631)</u>).

Acceptance criteria: The Sample solution is clear and not more intensely colored than alcohol or the Comparison solution.

• ACIDITY

Sample solution: To 2 mL of *Sample solution* prepared in the test for *Color of Solution*, add 3 mL of alcohol, 5 mL of carbon dioxide-free water, and 0.1 mL of bromocresol green TS.

Analysis: Titrate with 0.10 N sodium hydroxide.

Acceptance criteria: NMT 0.1 mL is required to produce a blue color.

ADDITIONAL REQUIREMENTS

• Packaging and Storage: Preserve in well-closed containers.

• USP REFERENCE STANDARDS (11)

USP Ethylparaben RS
USP Methylparaben RS

Auxiliary Information - Please check for your question in the FAQs before contacting USP.

| Topic/Question | Contact | Expert Committee |
|----------------|-------------------------------|--------------------------|
| ETHYLPARABEN | Documentary Standards Support | SE2020 Simple Excipients |

Chromatographic Database Information: Chromatographic Database

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. PF 36(5)

Current DocID: GUID-478339C6-47ED-41ED-AA0C-D0B5316E97FC_2_en-US

DOI: https://doi.org/10.31003/USPNF_M32170_02_01

DOI ref: n2a5x

