


Status: Currently Official on 15-Feb-2025
 Official Date: Official as of 01-Dec-2021
 Document Type: USP Monographs
 DocId: GUID-B155CCA7-9A30-44D0-AB7E-B0DB58F87B39_7_en-US
 DOI: https://doi.org/10.31003/USPNF_M41920_07_01
 DOI Ref: aw3of

© 2025 USPC
 Do not distribute

Iohexol



$C_{19}H_{26}I_3N_3O_9$ 821.14

1,3-Benzenedicarboxamide, 5-[acetyl(2,3-dihydroxypropyl)amino]-*N,N'*-bis(2,3-dihydroxypropyl)-2,4,6-triiodo;

N,N'-Bis(2,3-dihydroxypropyl)-5-[*N*-(2,3-dihydroxypropyl)acetamido]-2,4,6-triiodoisophthalamide CAS RN®: 66108-95-0; UNII: 4419T9MX03.

DEFINITION

Iohexol contains NLT 98.0% and NMT 102.0% of iohexol ($C_{19}H_{26}I_3N_3O_9$), calculated on the anhydrous basis.

IDENTIFICATION

- **A. SPECTROSCOPIC IDENTIFICATION TESTS (197), Infrared Spectroscopy:** 197K
- **B.** The retention times of the two principal peaks of the *Sample solution* correspond to those of the *System suitability solution*, as obtained in the test for *Organic Impurities*.

ASSAY

PROCEDURE

Sample: 500 mg of Iohexol

Sample solution: Transfer the *Sample* to a glass-stoppered, 125-mL conical flask. Add 25 mL of 1.25 N [sodium hydroxide](#) and 500 mg of powdered [zinc](#). Connect the flask to a reflux condenser, and reflux for 1 h. Cool the flask to room temperature, rinse the condenser with 20 mL of [water](#), disconnect the flask from the condenser, and pass the mixture through a filter. Rinse the flask and the filter thoroughly with small portions of [water](#), adding the rinsings to the filtrate. Add 5 mL of [glacial acetic acid](#).

Titrimetric system

(See [Titrimetry \(541\)](#).)

Mode: Direct titration

Titrant: [0.1 N silver nitrate VS](#)

Endpoint detection: Potentiometric

Analysis: Titrate the *Sample solution* with *Titrant*.

Calculate the percentage of iohexol ($C_{19}H_{26}I_3N_3O_9$) in the portion of Iohexol taken:

$$\text{Result} = [(V \times N \times F)/W] \times 100$$

V = *Titrant* volume consumed by the *Sample* (mL)

N = *Titrant* normality (mEq/mL)

F = equivalent weight of iohexol, 273.7 mg/mEq

W = *Sample* weight (mg)

Acceptance criteria: 98.0%–102.0% on the anhydrous basis

IMPURITIES

LIMIT OF IONIC COMPOUNDS

Rinse all glassware five times with distilled water.

Standard solution: 0.002 mg/mL of sodium chloride in [water](#)

Sample solution: 1 g of Iohexol in 50 mL of [water](#)

Analysis

Samples: *Standard solution* and *Sample solution*

Acceptance criteria: The specific conductance of the *Sample solution* is NMT that of the *Standard solution* (equivalent to 0.01% ionic compounds as sodium chloride).

Change to read:

• LIMIT OF FREE IODIDE

Sample: 5 g of Iohexol

Sample solution: Dissolve the *Sample* in 20 mL of [water](#).

Titrimetric system

(See [Titrimetry \(541\)](#).)

Mode: Direct titration

Titrant: 0.001 N [silver nitrate](#) ▲ (USP 1-Dec-2021)

Endpoint detection: Potentiometric

Analysis: Calculate the percentage of free iodide in the portion of Iohexol taken:

$$\text{Result} = [(V \times N \times F)/W] \times 100$$

V = *Titrant* volume consumed by the *Sample* (mL)

N = *Titrant* normality (mEq/mL)

F = equivalent weight of iodide, 126.9 mg/mEq

W = *Sample* weight (mg)

Acceptance criteria: NMT 0.001%

• ORGANIC IMPURITIES

Solution A: [Acetonitrile](#)

Solution B: [Water](#)

Mobile phase: See [Table 1](#).

Table 1

Time (min)	Solution A (%)	Solution B (%)
0	1	99
60	13	87

System suitability solution: 1.5 mg/mL of [USP Iohexol RS](#) and 0.0075 mg/mL of [USP Iohexol Related Compound A RS](#) in [water](#)

Sample solution: 1.5 mg/mL of Iohexol

Chromatographic system

(See [Chromatography \(621\)](#), [System Suitability](#).)

Mode: LC

Detector: UV 254 nm

Column: 4.6-mm × 25-cm; 5-μm packing [L1](#)

Flow rate: 1 mL/min

Injection volume: 10 μL

System suitability

Sample: *System suitability solution*

[NOTE—Iohexol may give two nonresolved peaks due to exo–endo isomerism. In addition, a small peak due to Iohexol usually appears at the leading edge of the first principal peak. This small peak has a retention time about 1.2 min less than the first principal peak. The relative retention times for the Iohexol related compound A, Iohexol endo-isomer, Iohexol exo-isomer, and O-alkylated compounds peaks are 0.85, 0.96, 1.0, and 1.1–1.4, respectively.]

Suitability requirements

Resolution: NLT 5.0 between iohexol related compound A and the exo-isomer (the second and greater peak) of iohexol

Analysis

Sample: *Sample solution*

Calculate the percentage of O-alkylated compounds and any other individual impurity in the portion of Iohexol taken. Exclude peaks with a relative retention time between 0.84 [relative to the endo-isomer of iohexol (first main peak)] and that of the endo-isomer of iohexol. Disregard any peak less than or equal to 0.03% of the principal peaks.

$$\text{Result} = (r_U/r_T) \times 100$$

r_U = peak response of each impurity

r_T = sum of all of the peak responses

Acceptance criteria

O-alkylated compounds: NMT 0.6%

Any individual impurity: NMT 0.1%

Total impurities excluding O-alkylated compounds: NMT 0.3%

Change to read:

• LIMIT OF 2-METHOXYETHANOL

Internal standard solution: ▲10 µg/mL▲ (USP 1-Dec-2021) of [secondary butyl alcohol](#) in [water](#)

▲**System suitability solution:** 5 µg/mL of [methanol](#) and 10 µg/mL each of [isopropyl alcohol](#) and [2-methoxyethanol](#) in *Internal standard solution*▲ (USP 1-Dec-2021)

Standard stock solution: ▲10 µg/mL of▲ (USP 1-Dec-2021) [2-methoxyethanol](#) in *Internal standard solution*

Standard solution: Transfer about 0.25 g of [USP Iohexol RS](#) and 1.0 mL of *Standard stock solution* to a headspace vial, and seal the vial with a septum and crimp cap.

Sample solution: Transfer about 0.25 g of Iohexol and 1.0 mL of *Internal standard solution* to a headspace vial, and seal the vial with a septum and crimp cap.

▲**Blank solution:** Transfer about 0.25 g of [USP Iohexol RS](#) and 1.0 mL of *Internal standard solution* to a headspace vial, and seal the vial with a septum and crimp cap.▲ (USP 1-Dec-2021)

Chromatographic system

(See [Chromatography \(621\)](#), [System Suitability](#).)

Mode: GC with suitable headspace autosampler

Detector: Flame ionization

Column: 0.53-mm × 30-m fused-silica; coated with a 1-µm phase [G16](#)

Temperatures

Autosampler: 105°

Needle: 130°–140°

Injection port: 150°

Detector: 200°

Column: See [Table 2](#).

Table 2

Initial Temperature (°)	Temperature Ramp (°/min)	Final Temperature (°)	Hold Time at Final Temperature (min)
40	—	40	3
40	8	100	1

Carrier gas: Helium

Flow rate: 11 mL/min

Injection volume: 1 mL of the headspace

System suitability

Sample: ▲*System suitability solution*▲ (USP 1-Dec-2021)

[NOTE—The typical relative retention times for methanol, isopropyl alcohol, secondary butyl alcohol, and 2-methoxyethanol are 0.5, 0.6, 1.0, and 1.9, respectively.]

Suitability requirements

Resolution: NLT 1.0 between methanol and isopropyl alcohol

Relative standard deviation: NMT 10.0% for the ratio of 2-methoxyethanol to ▲secondary butyl alcohol▲ (USP 1-Dec-2021)

Analysis

Samples: *Standard solution*, *Sample solution*, and ▲*Blank solution*▲ (USP 1-Dec-2021)

Calculate the amount of 2-methoxyethanol in the portion of Iohexol taken:

$$\text{Result} = [R_U / (R_S \text{▲} - R_B \text{▲})] \text{▲} \text{ (USP 1-Dec-2021)} \times (C_S / C_U)$$

R_U = peak response ratio of 2-methoxyethanol to ▲secondary butyl alcohol▲ (USP 1-Dec-2021) from the *Sample solution*

R_S = peak response ratio of 2-methoxyethanol to ▲secondary butyl alcohol▲ (USP 1-Dec-2021) from the *Standard solution*

▲ R_B = peak response ratio of 2-methoxyethanol to secondary butyl alcohol from the *Blank solution*▲ (USP 1-Dec-2021)

C_S = concentration of 2-methoxyethanol in the *Standard solution* (µg/mL)

C_U = concentration of Iohexol in the *Sample solution* (g/mL)

Acceptance criteria: NMT 20 µg/g of 2-methoxyethanol

Change to read:

• LIMIT OF 3-CHLOROPROPANE-1,2-DIOL

Standard solution: 0.025 mg/mL of [3-chloropropane-1,2-diol](#) in [ethyl acetate](#)

Sample solution: Transfer 1 g of Iohexol to a separator. Dissolve in 1 mL of [water](#). Extract 4 times with 2 mL of ethyl acetate, and combine the extracts. Dry the combined extracts with [anhydrous sodium sulfate](#). Filter, and wash the filter with a small amount of ethyl acetate. Combine the washings with the filtrate, and concentrate to a volume of 0.7 mL, using a warm water bath and a stream of nitrogen. Dilute with ethyl acetate to ▲1▲ (USP 1-Dec-2021) mL.

Chromatographic system

(See [Chromatography \(621\)](#), [System Suitability](#).)

Mode: GC

Detector: Flame ionization

Column: 0.32-mm × 30-m fused-silica capillary; bonded with a 1-µm layer of phase [G46](#)

Temperatures

Injection port: 230°

Detector: 250°

Column: See [Table 3](#).

Table 3

Initial Temperature (°)	Temperature Ramp (°/min)	Final Temperature (°)	Hold Time at Final Temperature (min)
80	—	80	2
80	15	275	—
275	—	275	2

Carrier gas: Helium

Flow rate: 1 mL/min

Injection volume: 2 µL

System suitability

Sample: *Standard solution*

[NOTE—The retention time of the 3-chloropropane-1,2-diol peak is about 8 min.]

Suitability requirements

Relative standard deviation: NMT 10.0%

Analysis

Samples: *Standard solution* and *Sample solution*

Acceptance criteria: The area of the principal peak from the *Sample solution* is NMT the area of the principal peak from the *Standard solution* (NMT 0.0025%).

• **LIMIT OF FREE AROMATIC AMINE**

Solution A: 3 mg/mL of [N-\(1-naphthyl\)ethylenediamine dihydrochloride](#) in a mixture of [propylene glycol](#) and [water](#) (70:30)

Standard stock solution: 10 µg/mL of [USP Iohexol Related Compound B RS](#) in [water](#)

Standard solution: Transfer 5 mL of [water](#) and 10.0 mL of the *Standard stock solution* to a 25-mL volumetric flask.

Sample solution: Transfer 200 mg of Iohexol to a 25-mL volumetric flask, add 15 mL of [water](#), and mix to dissolve.

Blank: Add 15 mL of [water](#) to a 25-mL volumetric flask.

Instrumental conditions

Mode: Vis

Analytical wavelength: 495 nm

Cell: 5 cm

Analysis

Samples: *Standard solution*, *Sample solution*, and *Blank*

In conducting the following steps, keep the flasks in iced water and protected as much as possible from light until all of the reagents have been added.

Treat the *Samples* as follows. Place the flask in an ice bath for 5 min. Add 1.5 mL of 6 N hydrochloric acid, and mix by swirling. Add 1.0 mL of sodium nitrite solution (20 mg/mL), and allow to stand in the ice bath for 4 min. Remove the flask from the ice bath, add 1.0 mL of sulfamic acid solution (40 mg/mL), and swirl gently until gas evolution ceases. [CAUTION—Considerable pressure is produced.] Add 1.0 mL of *Solution A*, dilute with water to volume, and allow to stand for 5 min. Measure the absorbance of the *Standard solution* and *Sample solution* against the *Blank*.

Acceptance criteria: The absorbance of the *Sample solution* is NMT that of the *Standard solution* (NMT 0.05% of free aromatic amine).

SPECIFIC TESTS

• **COLOR OF SOLUTION**

Sample solution: 647.2 mg/mL

Blank: [Water](#)

Instrumental conditions

Mode: UV-Vis

Analytical wavelengths: 400, 420, and 450 nm

Cell: 1 cm

Analysis

Samples: *Sample solution* and *Blank*

Pass the *Sample solution* through a filter of 0.22-µm pore size.

Determine the absorbances of the *Sample solution* against the *Blank*.

Acceptance criteria: See [Table 4](#).

Table 4

Wavelength (nm)	NMT (AU)
400	0.180
420	0.030
450	0.015

• [WATER DETERMINATION \(921\), Method I](#): NMT 4.0%

ADDITIONAL REQUIREMENTS

- **PACKAGING AND STORAGE:** Preserve in well-closed, light-resistant containers. Store at room temperature.
- **USP REFERENCE STANDARDS (11).**

[USP Iohexol RS](#)

[USP Iohexol Related Compound A RS](#)

5-(Acetylamino)-N,N'-bis(2,3-dihydroxypropyl)-2,4,6-triiodo-1,3-benzenedicarboxamide.
 $C_{16}H_{20}I_3N_3O_7$ 747.06

[USP Iohexol Related Compound B RS](#)

5-Amino-N,N'-bis(2,3-dihydroxypropyl)-2,4,6-triiodo-1,3-benzenedicarboxamide.
 $C_{14}H_{18}I_3N_3O_6$ 705.02

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
IOHEXOL	Documentary Standards Support	SM42020 Small Molecules 4

Chromatographic Database Information: [Chromatographic Database](#)

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. 46(4)

Current DocID: GUID-B155CCA7-9A30-44D0-AB7E-B0DB58F87B39_7_en-US

DOI: https://doi.org/10.31003/USPNF_M41920_07_01

DOI ref: [aw3of](#)

OFFICIAL