

Status: Currently Official on 16-Feb-2025
Official Date: Official Prior to 2013
Document Type: Reagents
DocId: GUID-E13CDB28-4821-4ABF-B593-302F9F854E03_1_en-US
DOI: https://doi.org/10.31003/USPNF_R3005_01_01
DOI Ref: 2fer0

© 2025 USPC
Do not distribute

Methyl Green–Iodomercurate Paper

—Immerse thin strips of suitable filter paper in a 40 g per L solution of methyl green, and allow to air-dry. Immerse the strips for 1 hour in a solution containing 140 g per L of potassium iodide and 200 g per L of mercuric iodide. Wash with water until the washings are practically colorless, and allow to air-dry. Store protected from light, and use within 48 hours.

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
METHYL GREEN–IODOMERCURATE PAPER	Margareth R.C. Marques Principal Scientific Liaison	HDQ Headquarters

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. Information currently unavailable

Current DocId: [GUID-E13CDB28-4821-4ABF-B593-302F9F854E03_1_en-US](#)

DOI: https://doi.org/10.31003/USPNF_R3005_01_01

DOI ref: [2fer0](#)

OFFICIAL