

Status: Currently Official on 16-Feb-2025

Official Date: Official as of 01-Aug-2017

Document Type: USP Monographs

DocId: GUID-1BA8E5C7-3870-4D15-AFDD-21FB0270DCAA_1_en-US

DOI: https://doi.org/10.31003/USPNF_M56180_01_01

DOI Ref: nt1d2

© 2025 USPC

Do not distribute

Neomycin and Polymyxin B Sulfates and Bacitracin Zinc Ophthalmic Ointment

DEFINITION

Neomycin and Polymyxin B Sulfates and Bacitracin Zinc Ophthalmic Ointment contains the equivalent of NLT 90.0% and NMT 140.0% of the labeled amounts of neomycin, polymyxin B, and bacitracin.

IDENTIFICATION

- [A. THIN-LAYER CHROMATOGRAPHIC IDENTIFICATION TEST \(201BNP\)](#): Meets the requirements

ASSAY

- **NEOMYCIN**

(See [Antibiotics—Microbial Assays \(81\)](#).)

Sample solution: Shake a portion of Ophthalmic Ointment in a separator with 50 mL of [ether](#). Extract with four 20-mL portions of *Buffer B.3*. Combine the aqueous extracts, and dilute with *Buffer B.3* to a suitable volume.

Analysis: Proceed as directed in the chapter. Dilute the *Sample solution* with *Buffer B.3* to obtain a *Test Dilution* having a neomycin concentration that is nominally equivalent to the median level of the standard.

Acceptance criteria: 90.0%–140.0%

- **POLYMYXIN B**

(See [Antibiotics—Microbial Assays \(81\)](#).)

Sample solution: Shake a portion of Ophthalmic Ointment with 50 mL of [ether](#) in a separator. Extract with four 25-mL portions of *Buffer B.6*. Combine the aqueous extracts, and dilute with *Buffer B.6* to a suitable volume.

Analysis: Proceed as directed in the chapter. Dilute the *Sample solution* with *Buffer B.6* to obtain a *Test Dilution* having a concentration that is nominally equivalent to the median level of the standard (10 Polymyxin B Units/mL). Add to each *Test Dilution* of the standard, a quantity of [USP Neomycin Sulfate RS](#), dissolved in *Buffer B.6*, to obtain the same concentration of neomycin as in the *Test Dilution* of the sample.

Acceptance criteria: 90.0%–140.0%

- **BACITRACIN**

(See [Antibiotics—Microbial Assays \(81\)](#).)

Sample solution: Shake a portion of Ophthalmic Ointment with 50 mL of [ether](#) in a separator. Extract with four 25-mL portions of [0.01 N hydrochloric acid](#). Combine the acid extracts, and dilute with [0.01 N hydrochloric acid](#) to a suitable volume.

Analysis: Proceed as directed in the chapter. Dilute the *Sample solution* with *Buffer B.1* to obtain a *Test Dilution* having a concentration that is nominally equivalent to the median level of the standard (1.0 Bacitracin Unit/mL). If the *Sample solution* has a concentration of less than 100 Bacitracin Units/mL, add hydrochloric acid to each *Test Dilution* of the standard to obtain the same concentration of hydrochloric acid as in the *Test Dilution* of the sample.

Acceptance criteria: 90.0%–140.0%

SPECIFIC TESTS

- [STERILITY TESTS \(71\)](#): Meets the requirements
- **OTHER REQUIREMENTS:** It meets the requirements for Particulate and Foreign Matter and Container Contents in [Ophthalmic Products—Quality Tests \(771\)](#), [Drug Product Quality, Universal Tests, Particulate and Foreign Matter](#) and [Container Contents](#).

ADDITIONAL REQUIREMENTS

- **PACKAGING AND STORAGE:** Preserve in collapsible ophthalmic ointment tubes. Store at controlled room temperature.

- [USP REFERENCE STANDARDS \(11\)](#):

[USP Bacitracin Zinc RS](#)

[USP Neomycin Sulfate RS](#)

[USP Polymyxin B Sulfate RS](#)

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
NEOMYCIN AND POLYMYXIN B SULFATES AND BACITRACIN ZINC OPHTHALMIC OINTMENT	Julie Zhang Associate Science & Standards Liaison	BIO42020 Biologics Monographs 4 - Antibiotics
REFERENCE STANDARD SUPPORT	RS Technical Services RSTECH@usp.org	BIO42020 Biologics Monographs 4 - Antibiotics

Chromatographic Database Information: [Chromatographic Database](#)

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. PF 39(5)

Current DocID: GUID-1BA8E5C7-3870-4D15-AFDD-21FB0270DCAA_1_en-US

DOI: https://doi.org/10.31003/USPNF_M56180_01_01

DOI ref: [nt1d2](#)

OFFICIAL