


Status: Currently Official on 16-Feb-2025  
Official Date: Official Prior to 2013  
Document Type: USP Monographs  
DocId: GUID-FDC010A2-F2C8-4431-AF61-D37522E5AB74\_1\_en-US  
DOI: [https://doi.org/10.31003/USPNF\\_M73550\\_01\\_01](https://doi.org/10.31003/USPNF_M73550_01_01)  
DOI Ref: 457j2

© 2025 USPC  
Do not distribute

## Riboflavin 5'-Phosphate Sodium



$C_{17}H_{20}N_4NaO_9P \cdot 2H_2O$  514.36

$C_{17}H_{20}N_4NaO_9P$  478.33

Riboflavin 5'-(dihydrogen phosphate), monosodium salt, dihydrate;

Riboflavin 5'-(sodium hydrogen phosphate), dihydrate

Anhydrous CAS RN®: 130-40-5; UNII: 957E53WV42.

### DEFINITION

Riboflavin 5'-Phosphate Sodium contains NLT the equivalent of 73.0% and NMT the equivalent of 79.0% of riboflavin ( $C_{17}H_{20}N_4O_6$ ), calculated on the dried basis.

### IDENTIFICATION

#### • A. COLOR AND FLUORESCENCE OF SOLUTION

**Sample solution:** 0.01 mg/mL in water

**Analysis:** Alternately expose to transmitted light and long-wavelength UV light.

**Acceptance criteria:** The solution is pale greenish yellow by transmitted light, and it exhibits an intense yellowish green fluorescence by reflected light under long-wavelength UV light that disappears upon the addition of mineral acids or alkalies.

#### • B. IDENTIFICATION TESTS—GENERAL, Sodium(191)andPhosphate(191)

**Sample solution:** To 0.5 g add 10 mL of nitric acid. Evaporate the mixture on a water bath to dryness, and ignite the residue until the carbon is removed. Dissolve the residue in 5 mL of water, and filter.

**Acceptance criteria:** The *Sample solution* meets the requirements.

### ASSAY

#### • PROCEDURE

[**NOTE**—Conduct the assay so that all solutions are protected from actinic light at all stages, preferably by using low-actinic glassware.]

**Standard solution:** 0.35  $\mu$ g/mL prepared as follows. Transfer 35 mg of [USP Riboflavin RS](#) to a 250-mL conical flask. Add 20 mL of pyridine and 75 mL of water, and dissolve the riboflavin by frequent shaking. Transfer the solution to a 1000-mL volumetric flask, and dilute with water to volume. Transfer 10.0 mL of this solution to a second 1000-mL volumetric flask, add sufficient 0.1 N sulfuric acid (about 4 mL) so that the final pH of the solution is between 5.9 and 6.1, and dilute with water to volume.

**Sample solution:** Transfer 50 mg of Riboflavin 5'-Phosphate Sodium to a 250-mL conical flask. Add 20 mL of pyridine and 75 mL of water, and dissolve the riboflavin by frequent shaking. Transfer the solution to a 1000-mL volumetric flask, and dilute with water to volume. Transfer 10.0 mL of this solution to a second 1000-mL volumetric flask, add sufficient 0.1 N sulfuric acid (about 4 mL) so that the final pH of the solution is between 5.9 and 6.1, and dilute with water to volume.

**Blank:** Proceed as directed in the *Analysis*, omitting the test specimen.

### Instrumental conditions

(See [Fluorescence Spectroscopy \(853\)](#).)

**Mode:** Fluorescence

**Excitation wavelength:** 440 nm

**Emission wavelength:** 530 nm

**Analysis****Samples:** Standard solution, Sample solution, and Blank

Determine the maximum fluorescence intensities of the solutions against the Blank.

Calculate the percentage of riboflavin ( $C_{17}H_{20}N_4O_6$ ) in the portion of Riboflavin 5'-Phosphate Sodium taken:

$$\text{Result} = (I_U/I_S) \times (C_S/C_U) \times 100$$

$I_U$  = fluorescence intensity from the Sample solution

$I_S$  = fluorescence intensity from the Standard solution

$C_S$  = concentration of [USP Riboflavin RS](#) in the Standard solution ( $\mu\text{g/mL}$ )

$C_U$  = concentration of Riboflavin 5'-Phosphate Sodium in the Sample solution ( $\mu\text{g/mL}$ )

**Acceptance criteria:** 73.0%–79.0% on the dried basis

**IMPURITIES****• FREE PHOSPHATE**

**Acid molybdate solution:** Prepare a 70-mg/mL solution of ammonium molybdate in water. Dilute 25 mL of this solution with water to 200 mL. To this dilution slowly add 25 mL of 7.5 N sulfuric acid.

**Ferrous sulfate solution:** 100 mg/mL of ferrous sulfate in 0.15 N sulfuric acid, prepared just before use

**Standard solution:** 44.0  $\mu\text{g/mL}$  of monobasic potassium phosphate in water.

**Sample solution:** 3 mg/mL of Riboflavin 5'-Phosphate Sodium in water

**Blank:** Water

**Instrumental conditions**

(See [Ultraviolet-Visible Spectroscopy \(857\)](#).)

**Mode:** UV-Vis

**Analytical wavelength:** 700 nm

**Cell:** 1 cm

**Analysis****Samples:** Standard solution, Sample solution, and Blank

Transfer 10.0 mL each of the Standard solution, Sample solution, and Blank to separate 50-mL conical flasks. Add 10.0 mL of Acid molybdate solution and 5.0 mL of Ferrous sulfate solution to each flask, and mix. Determine the absorbances of the solutions against that of the Blank.

**Acceptance criteria:** NMT 1% as  $\text{PO}_4$ . The absorbance of the Sample solution is NMT that of the Standard solution.

**• FREE RIBOFLAVIN AND RIBOFLAVIN DIPHOSPHATES**

[NOTE—Conduct this test so that all solutions are protected from actinic light at all stages, preferably by using low-actinic glassware.]

**Mobile phase:** Methanol and 0.054 M monobasic potassium phosphate (15:85)

**System suitability solution:** Prepare a 2-mg/mL solution of [USP Phosphated Riboflavin RS](#) in water. Dilute this solution with Mobile phase to 160  $\mu\text{g/mL}$ .

**Standard solution:** Transfer 60 mg of [USP Riboflavin RS](#) to a 250-mL volumetric flask. Dissolve carefully in 1 mL of hydrochloric acid, and dilute with water to volume. Dilute with Mobile phase to 9.6  $\mu\text{g/mL}$ .

**Sample solution:** Prepare a 2-mg/mL solution of Riboflavin 5'-Phosphate Sodium in water. Dilute this solution with Mobile phase to 160  $\mu\text{g/mL}$ .

**Chromatographic system**

(See [Chromatography \(621\), System Suitability](#).)

**Mode:** LC

**Detector:** Fluorometer

**Excitation wavelength:** 440 nm

**Emission wavelength:** 530 nm (monochromator-based detector) or 470 nm (filtered-type detector)

**Column:** 3.9-mm  $\times$  30-cm; packing L1

**Flow rate:** 2 mL/min

**Injection size:** 100  $\mu\text{L}$

**System suitability**

**Sample:** System suitability solution

The retention time of riboflavin 5'-monophosphate is 20–25 min. The approximate relative retention times for the components are listed in [Table 1](#).

**Table 1**

Name	Relative Retention Time
Riboflavin 3'4'-diphosphate	0.23
Riboflavin 3'5'-diphosphate	0.39
Riboflavin 4'5'-diphosphate	0.58
Riboflavin 3'-monophosphate	0.70
Riboflavin 4'-monophosphate	0.87
Riboflavin 5'-monophosphate	1.00
Riboflavin	1.63

**Suitability requirements**

**Resolution:** NLT 1.0 between riboflavin 4'-monophosphate and riboflavin 5'-monophosphate

**Relative standard deviation:** NMT 1.5% for riboflavin 5'-monophosphate

**Analysis**

**Samples:** Standard solution and Sample solution

Measure the peak areas of the solutions. Identify the peaks to be measured in the *Sample solution* by comparison with the chromatogram obtained from the *System suitability solution*.

Calculate the percentage of free riboflavin in the portion of Riboflavin 5'-Phosphate Sodium taken:

$$\text{Result} = (r_u/r_s) \times (C_s/C_u) \times 100$$

$r_u$  = peak area of free riboflavin from the *Sample solution*

$r_s$  = peak area of riboflavin from the *Standard solution*

$C_s$  = concentration of [USP Riboflavin RS](#) in the *Standard solution* ( $\mu\text{g/mL}$ )

$C_u$  = concentration of Riboflavin 5'-Phosphate Sodium in the *Sample solution* ( $\mu\text{g/mL}$ )

Calculate the percentage of riboflavin diphosphates in the portion of Riboflavin 5'-Phosphate Sodium taken:

$$\text{Result} = (r_u/r_s) \times (C_s/C_u) \times 100$$

$r_u$  = sum of the peak areas for any of the 3 riboflavin diphosphates from the *Sample solution*

$r_s$  = peak area of riboflavin from the *Standard solution*

$C_s$  = concentration of [USP Riboflavin RS](#) in the *Standard solution* ( $\mu\text{g/mL}$ )

$C_u$  = concentration of Riboflavin 5'-Phosphate Sodium in the *Sample solution* ( $\mu\text{g/mL}$ )

**Acceptance criteria**

**Free riboflavin:** NMT 6.0% on the dried basis

**Riboflavin diphosphates, as riboflavin:** NMT 6.0% on the dried basis

**• LIMIT OF LUMIFLAVIN**

**Alcohol-free chloroform:** Shake 20 mL of chloroform gently but thoroughly with 20 mL of water for 3 min, draw off the chloroform layer, and wash twice more with 20-mL portions of water. Finally, pass the chloroform through a dry filter paper, shake it for 5 min with 5 g of powdered anhydrous sodium sulfate, allow the mixture to stand for 2 h, and decant or filter the clear chloroform.

**Sample solution:** Shake 35 mg of Riboflavin 5'-Phosphate Sodium with 10 mL of *Alcohol-free chloroform* for 5 min, and filter.

**Blank:** Alcohol-free chloroform

#### Instrumental conditions

(See [Ultraviolet-Visible Spectroscopy \(857\)](#).)

**Analytical wavelength:** 440 nm

**Cell:** 1 cm

#### Analysis

**Samples:** *Sample solution* and *Blank*

Measure the absorbances of the *Sample solution* and *Blank*. Correct the absorbance of the *Sample solution* with the *Blank*.

**Acceptance criteria:** The absorbance is NMT 0.025.

#### SPECIFIC TESTS

- [RESIDUE ON IGNITION \(281\)](#): NMT 25.0%
- [OPTICAL ROTATION, Specific Rotation \(781S\)](#).

**Sample solution:** 15 mg/mL in 5 N hydrochloric acid. Use the solution within 15 min.

**Acceptance criteria:** +37.0° to +42.0°

- [pH \(791\)](#).

**Sample solution:** 10 mg/mL solution

**Acceptance criteria:** 5.0–6.5

- [Loss on Drying \(731\)](#): Dry a sample under vacuum over phosphorus pentoxide at 100° for 5 h; it loses NMT 7.5% of its weight.

#### ADDITIONAL REQUIREMENTS

- **PACKAGING AND STORAGE:** Preserve in tight, light-resistant containers.

- [USP REFERENCE STANDARDS \(11\)](#).

[USP Phosphated Riboflavin RS](#)

[USP Riboflavin RS](#)

**Auxiliary Information** - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
RIBOFLAVIN 5'-PHOSPHATE SODIUM	<a href="#">Natalia Davydova</a> Scientific Liaison	NBDS2020 Non-botanical Dietary Supplements
REFERENCE STANDARD SUPPORT	RS Technical Services <a href="mailto:RSTECH@usp.org">RSTECH@usp.org</a>	NBDS2020 Non-botanical Dietary Supplements

**Chromatographic Database Information:** [Chromatographic Database](#)

#### Most Recently Appeared In:

Pharmacopeial Forum: Volume No. Information currently unavailable

**Current DocID: GUID-FDC010A2-F2C8-4431-AF61-D37522E5AB74\_1\_en-US**

**DOI:** [https://doi.org/10.31003/USPNF\\_M73550\\_01\\_01](https://doi.org/10.31003/USPNF_M73550_01_01)

**DOI ref:** [457j2](#)