

Status: Currently Official on 16-Feb-2025
 Official Date: Official as of 01-Aug-2023
 Document Type: Reagents
 DocId: GUID-ED3B9FD9-AED0-4703-BEFE-CD6130BA7C98_2_en-US
 DOI: https://doi.org/10.31003/USPNF_R8061_02_01
 DOI Ref: g4dec

© 2025 USPC
 Do not distribute

Add the following:

0.01 N Silver Nitrate VS

Transfer 100 mL of 0.1 N silver nitrate VS to a 1000-mL volumetric flask. Dilute with water to volume. Otherwise, dissolve about 1.75 g of silver nitrate in 1000 mL of water.

Standardization

See [Volumetric Solutions, 1. Introduction](#).

See [Titrimetry \(541\)](#).

Standardize by one of the following procedures. [NOTE—Other standardization procedures may be used. See [Volumetric Solutions, 2. Preparation and Standardization, 2.3 Standardization](#).]

Standardization with visual endpoint: Transfer 0.50 mL of [0.1 N sodium chloride VS](#) to a 150-mL beaker containing 50 mL of [water](#), and add 50 mL of [alcohol](#), 5 mL of [acetic acid](#), and 0.5 mL of [eosin Y TS](#). Stir, preferably with a magnetic stirrer, and titrate with the silver nitrate solution.

Standardization with potentiometric endpoint: Transfer 0.50 mL of [0.1 N sodium chloride VS](#) to a 150-mL beaker containing 50 mL of [water](#), and add 50 mL of [alcohol](#), 5 mL of [acetic acid](#). Stir, preferably with a magnetic stirrer, and titrate with the silver nitrate solution using a combined silver electrode.

[NOTE—If this volumetric solution is used in a qualitative application such as pH adjustment, dissolution medium, or diluent, its standardization is not required.]▲ (USP 1-Aug-2023)

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
0.01 N SILVER NITRATE VS	Margareth R.C. Marques Principal Scientific Liaison	HDQ Headquarters

Most Recently Appeared In:
 Pharmacopeial Forum: Volume No. 48(2)

Current DocID: [GUID-ED3B9FD9-AED0-4703-BEFE-CD6130BA7C98_2_en-US](#)

DOI: https://doi.org/10.31003/USPNF_R8061_02_01

DOI ref: [g4dec](#)