Status: Currently Official on 14-Feb-2025
Official Date: Official as of 01-May-2018
Document Type: USP Monographs
DocId: GUID-38FF7566-60C8-44B2-B602-86BF79518198_3_en-US
DOI: https://doi.org/10.31003/USPNF_M15140_03_01
DOI Ref: 52kt3

© 2025 USPC Do not distribute

Chloramphenicol Injection

» Chloramphenicol Injection is a sterile solution of Chloramphenicol in one or more suitable solvents. It contains not less than 90.0 percent and not more than 115.0 percent of the labeled amount of $C_{11}H_{12}Cl_2N_2O_5$. It may contain suitable buffers.

Packaging and storage—Preserve in single-dose or in multiple-dose containers.

Labeling—Label it to indicate that it is for veterinary use only.

USP REFERENCE STANDARDS (11)-

USP Chloramphenicol RS

Identification—The retention time of the major peak in the chromatogram of the *Assay preparation* corresponds to that in the chromatogram of the *Standard preparation* as obtained in the *Assay*.

BACTERIAL ENDOTOXINS TEST (85) —It contains not more than 0.2 USP Endotoxin Unit per mg of chloramphenicol.

<u>STERILITY TESTS (71)</u>—It meets the requirements when tested as directed for *Membrane Filtration* under *Test for Sterility of the Product to be Examined*, 1 mL from each container being transferred directly to the membrane filter.

PH (791): between 5.0 and 8.0, in a solution diluted with water (1:1).

Other requirements—It meets the requirements under <u>Injections and Implanted Drug Products (1)</u>.

Assay-

Mobile phase, Standard preparation, and Chromatographic system—Prepare as directed in the Assay under Chloramphenicol.

Assay preparation—Transfer an accurately measured volume of Chloramphenicol Injection, equivalent to about 200 mg of chloramphenicol, to a 100-mL volumetric flask, add *Mobile phase* to volume, and mix. Transfer 4.0 mL of the resulting solution to a 100-mL volumetric flask, dilute with *Mobile phase* to volume, and mix. Filter this solution through a 0.5-µm or finer porosity filter.

Procedure—Proceed as directed for <u>Procedure</u> in the <u>Assay</u> under <u>Chloramphenicol</u>. Calculate the quantity, in mg, of C₁₁H₁₂Cl₂N₂O₅ in each mL of the Injection taken by the formula:

 $2.5(C/V)(r_{11}/r_{s})$

in which V is the volume, in mL, of Injection taken, and the other terms are as defined therein.

Auxiliary Information - Please check for your question in the FAQs before contacting USP.

Topic/Question	Contact	Expert Committee
CHLORAMPHENICOL INJECTION	Documentary Standards Support	SM32020 Small Molecules 3

Chromatographic Database Information: Chromatographic Database

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. Information currently unavailable

Current DocID: GUID-38FF7566-60C8-44B2-B602-86BF79518198_3_en-US Previous DocID: GUID-38FF7566-60C8-44B2-B602-86BF79518198_1_en-US

DOI: https://doi.org/10.31003/USPNF_M15140_03_01

DOI ref: 52kt3