

Status: Currently Official on 14-Feb-2025  
Official Date: Official Prior to 2013  
Document Type: Reagents  
DocId: GUID-2A2ADBC3-B70F-4F74-ABD1-2530E2B3BA1A\_1\_en-US  
DOI: [https://doi.org/10.31003/USPNF\\_R3159\\_01\\_01](https://doi.org/10.31003/USPNF_R3159_01_01)  
DOI Ref: f7p0m

© 2025 USPC  
Do not distribute

## Dinitrophenylhydrazine TS

—Carefully mix 10 mL of water and 10 mL of sulfuric acid, and cool. To the mixture, contained in a glass-stoppered flask, add 2 g of 2,4-dinitrophenylhydrazine, and shake until dissolved. To the solution add 35 mL of water, mix, cool, and filter.

**Auxiliary Information** - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
DINITROPHENYLHYDRAZINE TS	<a href="#">Margareth R.C. Marques</a> Principal Scientific Liaison	HDQ Headquarters

**Most Recently Appeared In:**

Pharmacopeial Forum: Volume No. Information currently unavailable

**Current DocID:** [GUID-2A2ADBC3-B70F-4F74-ABD1-2530E2B3BA1A\\_1\\_en-US](#)

**DOI:** [https://doi.org/10.31003/USPNF\\_R3159\\_01\\_01](https://doi.org/10.31003/USPNF_R3159_01_01)

**DOI ref:** [f7p0m](#)

OFFICIAL