

Status: Currently Official on 15-Feb-2025
Official Date: Official Prior to 2013
Document Type: Reagents
DocId: GUID-7DD1281C-49CB-416F-8415-E9A21E4DCF6A_1_en-US
DOI: https://doi.org/10.31003/USPNF_R3177_01_01
DOI Ref: tk53m

© 2025 USPC
Do not distribute

Gelatin TS

(for the assay of *Corticotropin Injection*)—Dissolve 340 g of acid-treated precursor gelatin (Type A) in water to make 1000 mL. Heat the solution in an autoclave at 115° for 30 minutes after the exhaust line temperature has reached 115°. Cool the solution, and add 10 g of phenol and 1000 mL of water. Store in tight containers in a refrigerator.

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
GELATIN TS	Margareth R.C. Marques Principal Scientific Liaison	HDQ Headquarters

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. Information currently unavailable

Current DocId: [GUID-7DD1281C-49CB-416F-8415-E9A21E4DCF6A_1_en-US](#)

DOI: https://doi.org/10.31003/USPNF_R3177_01_01

DOI ref: [tk53m](#)

OFFICIAL