

Status: Currently Official on 15-Feb-2025  
Official Date: Official Prior to 2013  
Document Type: Reagents  
DocId: GUID-D0D2F8C7-C3B2-4CFF-B53C-56841A0C897D\_1\_en-US  
DOI: [https://doi.org/10.31003/USPNF\\_R3558\\_01\\_01](https://doi.org/10.31003/USPNF_R3558_01_01)  
DOI Ref: bzg7i

© 2025 USPC  
Do not distribute

## L33

—Packing having the capacity to separate dextrans by molecular size over a range of 4,000 to 500,000 Da. It is spherical, silica-based, and processed to provide pH stability. [NOTE—Available as TSK-GEL G4000SWxl from Tosoh Bioscience ([www.tosohbioscience.com](http://www.tosohbioscience.com)).]

**Auxiliary Information** - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
L33	<a href="#">Margareth R.C. Marques</a> Principal Scientific Liaison	HDQ Headquarters

**Most Recently Appeared In:**

Pharmacopeial Forum: Volume No. Information currently unavailable

**Current DocID:** [GUID-D0D2F8C7-C3B2-4CFF-B53C-56841A0C897D\\_1\\_en-US](#)

**DOI:** [https://doi.org/10.31003/USPNF\\_R3558\\_01\\_01](https://doi.org/10.31003/USPNF_R3558_01_01)

**DOI ref:** [bzg7i](#)

