

Status: Currently Official on 15-Feb-2025  
Official Date: Official Prior to 2013  
Document Type: Reagents  
DocId: GUID-5C25B766-58EC-46E8-9CE2-DD42E9C97F46\_1\_en-US  
DOI: [https://doi.org/10.31003/USPNF\\_R3201\\_01\\_01](https://doi.org/10.31003/USPNF_R3201_01_01)  
DOI Ref: 4f7tn

© 2025 USPC  
Do not distribute

## Litmus TS

—Digest 25 g of powdered litmus with three successive 100-mL portions of boiling alcohol, continuing each extraction for about 1 hour. Filter, wash with alcohol, and discard the alcohol filtrate. Macerate the residue with about 25 mL of cold water for 4 hours, filter, and discard the filtrate. Finally digest the residue with 125 mL of boiling water for 1 hour, cool, and filter.

**Auxiliary Information** - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
LITMUS TS	<a href="#">Margareth R.C. Marques</a> Principal Scientific Liaison	HDQ Headquarters

**Most Recently Appeared In:**

Pharmacopeial Forum: Volume No. Information currently unavailable

**Current DocId:** [GUID-5C25B766-58EC-46E8-9CE2-DD42E9C97F46\\_1\\_en-US](#)

**DOI:** [https://doi.org/10.31003/USPNF\\_R3201\\_01\\_01](https://doi.org/10.31003/USPNF_R3201_01_01)

**DOI ref:** [4f7tn](#)

OFFICIAL