

Status: Currently Official on 16-Feb-2025  
Official Date: Official Prior to 2013  
Document Type: Reagents  
DocId: GUID-677B59A2-8533-43CD-A92F-A9A0DF230DB5\_1\_en-US  
DOI: [https://doi.org/10.31003/USPNF\\_R3233\\_01\\_01](https://doi.org/10.31003/USPNF_R3233_01_01)  
DOI Ref: k2y8k

© 2025 USPC  
Do not distribute

## Molybdo-phosphotungstate TS (*Folin-Denis Reagent*)—

To about 350 mL of water contained in a round-bottom flask add 50 g of sodium tungstate, 12 g of phosphomolybdic acid, and 25 mL of phosphoric acid. Boil the mixture under a reflux condenser for 2 hours, then cool, dilute with water to 500 mL, and mix. Store in tight containers, protected from light, and in a cold place.

**Auxiliary Information** - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
MOLYBDO-PHOSPHOTUNGSTATE TS (FOLIN-DENIS REAGENT)	<a href="#">Margareth R.C. Marques</a> Principal Scientific Liaison	HDQ Headquarters

**Most Recently Appeared In:**

Pharmacopeial Forum: Volume No. Information currently unavailable

**Current DocId:** [GUID-677B59A2-8533-43CD-A92F-A9A0DF230DB5\\_1\\_en-US](#)

**DOI:** [https://doi.org/10.31003/USPNF\\_R3233\\_01\\_01](https://doi.org/10.31003/USPNF_R3233_01_01)

**DOI ref:** [k2y8k](#)

