

Status: Currently Official on 14-Feb-2025
Official Date: Official Prior to 2013
Document Type: Reagents
DocId: GUID-C069FC1D-6E7A-4768-B50F-F90F195F91CA_1_en-US
DOI: https://doi.org/10.31003/USPNF_R3120_01_01
DOI Ref: tt9nj

© 2025 USPC
Do not distribute

p-Bromoaniline TS

—Add 8 g of p-bromoaniline to a mixture of 380 mL of thiourea-saturated glacial acetic acid, 10 mL of sodium chloride solution (1 in 5), 5 mL of oxalic acid solution (1 in 20), and 5 mL of dibasic sodium phosphate solution (1 in 10) in a low-actinic glass bottle. Mix, and allow to stand overnight before using. Protect from light, and use within 7 days.

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
P-BROMOANILINE TS	Margareth R.C. Marques Principal Scientific Liaison	HDQ Headquarters

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. Information currently unavailable

Current DocId: [GUID-C069FC1D-6E7A-4768-B50F-F90F195F91CA_1_en-US](#)

DOI: https://doi.org/10.31003/USPNF_R3120_01_01

DOI ref: [tt9nj](#)

OFFICIAL