

Status: Currently Official on 16-Feb-2025
Official Date: Official Prior to 2013
Document Type: USP Monographs
DocId: GUID-4C7D1E2E-7310-424B-9CA0-6CCED1C551E0_1_en-US
DOI: https://doi.org/10.31003/USPNF_M81743_01_01
DOI Ref: qug7d

© 2025 USPC
Do not distribute

Tetracycline Boluses

» Tetracycline Boluses contain the equivalent of not less than 90.0 percent and not more than 120.0 percent of the labeled amount of tetracycline hydrochloride ($C_{22}H_{24}N_2O_8 \cdot HCl$).

Packaging and storage—Preserve in tight containers.

Labeling—Label Boluses to indicate that they are intended for veterinary use only.

USP REFERENCE STANDARDS (11).—
[USP Tetracycline Hydrochloride RS](#)

Identification—Shake a suitable quantity of finely powdered Boluses with methanol to obtain a solution containing the equivalent of 1 mg of tetracycline hydrochloride per mL, and filter. Using the filtrate as the *Test Solution*, proceed as directed under [Identification—Tetracyclines \(193\)](#).

UNIFORMITY OF DOSAGE UNITS (905): meet the requirements for *Weight Variation*.

LOSS ON DRYING (731).—Dry about 100 mg of finely powdered Boluses, accurately weighed, in a capillary-stoppered bottle in vacuum at a pressure not exceeding 5 mm of mercury at 60° for 3 hours: it loses not more than 3.0% of its weight, or where the Boluses have a diameter of greater than 15 mm, it loses not more than 6.0% of its weight.

Assay—Transfer not less than 2 Boluses to a high-speed blender jar containing an accurately measured volume of 0.1 N hydrochloric acid, so that the solution so obtained contains not less than 150 µg of tetracycline hydrochloride per mL, and blend for about 3 to 5 minutes. Proceed as directed for tetracycline under [Antibiotics—Microbial Assays \(81\)](#), using an accurately measured volume of this solution diluted quantitatively and stepwise with water to obtain a *Test Dilution* having a concentration assumed to be equal to the median dose level of the Standard.

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
TETRACYCLINE BOLUSES	Julie Zhang Associate Science & Standards Liaison	BIO42020 Biologics Monographs 4 - Antibiotics
REFERENCE STANDARD SUPPORT	RS Technical Services RSTECH@usp.org	BIO42020 Biologics Monographs 4 - Antibiotics

Chromatographic Database Information: [Chromatographic Database](#)

Most Recently Appeared In:
Pharmacopeial Forum: Volume No. Information currently unavailable

Current DocID: GUID-4C7D1E2E-7310-424B-9CA0-6CCED1C551E0_1_en-US
DOI: https://doi.org/10.31003/USPNF_M81743_01_01
DOI ref: [qug7d](#)