

Status: Currently Official on 17-Feb-2025
Official Date: Official Prior to 2013
Document Type: Reagents
DocId: GUID-535F14AB-87EC-4E00-B666-8CC07B4B4C36_1_en-US
DOI: https://doi.org/10.31003/USPNF_R3348_01_01
DOI Ref: sl0mh

© 2025 USPC
Do not distribute

Zinc Uranyl Acetate TS

—Dissolve 50 g of uranyl acetate in a mixture of 15 mL of glacial acetic acid and water to make 500 mL. Then dissolve 150 g of zinc acetate in a mixture of 15 mL of glacial acetic acid and water to make 500 mL. Mix the two solutions, allow to stand overnight, and pass through a dry filter, if necessary.

Auxiliary Information - Please [check for your question in the FAQs](#) before contacting USP.

Topic/Question	Contact	Expert Committee
ZINC URANYL ACETATE TS	Margareth R.C. Marques Principal Scientific Liaison	HDQ Headquarters

Most Recently Appeared In:

Pharmacopeial Forum: Volume No. Information currently unavailable

Current DocId: [GUID-535F14AB-87EC-4E00-B666-8CC07B4B4C36_1_en-US](#)

DOI: https://doi.org/10.31003/USPNF_R3348_01_01

DOI ref: [sl0mh](#)

OFFICIAL